Spekdiktemeter Renco
Renco lean indicator
Renco Fettschichtmesser
Appareil pour mesurer l'épaisseur de lard Renco
Spessimetro per grasso
1603446

NL	Spekdiktemeter Renco	
	Gebruiksaanwijzing	2
EN	Renco lean indicator	
	Manual	5
DE	Renco Fettschichtmesser	
	Mode d'emploi	7
FR	Appareil pour mesurer l'épaisseur de lard Renco	
	Mode d'emploi	8
IT	Spessimetro per grasso	
	Istruzioni per l'uso	11
PL	Renco Lean	
	Podręcznik	13

The
Schippers
 Group

Let's beat AMR with HyCare

NL Klantenservice	+31 (0)497 339 787
BE Klantenservice/Service Clientèle	+32 (0)14-820713
DE Kundenservice	+49 (0)2833-923630
DK Kundeservice	+45 89884187
IT Servizio clienti	035-4490369
ES Servicio al cliente	+34 931816433
UK Customer Service	+44 (0)1733592049
CA Customer Service	+1 866 995-7771
FR Service Clientèle	+33 (0)2 99 61 40 40
EX Customer Service Export	+31 (0)497 700 278

Spekdiktemeter Renco

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
Renco lean indicator

Renco lean indicator

All hogs have a genetic potential for three layers of fat. As the hog grows, fat is deposited in the top two layers. The third layer will not have a fat deposit until later in the hog's life. When the third layer fills with fat, you can measure the dept of the three layers of fat. **Do not attempt to measure three layers of fat over the shoulder or the ham.** The resulting measurements may not be accurate.

1. How to use the Renco lean indicator

1. Use the Renco charger supplied to charge the lean indicator battery for about 12 hours before the first use. Charging the battery after each day's use for three times longer than the time it was used will insure long battery life. When the battery becomes discharged the display will either begin flashing, or remain off.

2. The most accurate spot in which to measure backfat can be found by walking your fingers forward along the hog's flank until you can feel the last rib, then place the probe on this site B, 65mm (2.6 inches) from either side of the backbone. Saturate the skin with a light oil, like scan gel or paraffin oil on the spot where apply the lean indicator. Insufficient amount of fluids is the most common cause for inaccurate readings. Good skin contact with the probe and skin must be made for the ultrasound to enter the hog and reflect back to the lean indicator.

To make a measurement, turn the instrument on by pressing the POWER push button. Apply the probe to the test site, and using a gentle pressing/twisting action, move the probe slightly to force out any air bubbles between the skin and probe. **It is very important that the probe be kept perpendicular to the back. Errors may result if probe is at an angle.** Many have found that probe is most easily managed if held between the index and middle fingers.

When the READ light, to the left of the main display, lights the reading is complete. If the READ light does not come on, then skin contact is poor. Use more oil. Make a greater effort to remove air bubbles. Be sure probe is flat against the skin and perpendicular to the animal's back. Consistency of probe placement is of great importance in obtaining comparative measurements.

3. Turn the instrument on. The red READ indicator will light when the lean indicator detects and measures, the total of the fat layers. The lean indicator attempts a new measurement several times each second. A reading will remain until a new, different reading is detected.

Small irregularities in the fat may sometimes block the beam of sound, causing the lean indicator to give low readings.

Moving the probe slightly will allow the beam to bypass, and the correct reading to appear. **(Always keep the probe perpendicular to the animal's back)**

Fat layers are not uniformly thick. Thus, averaging two or three readings at adjacent sites insures greatest accuracy. The initial reading that occurs when the POWER button is pressed has no meaning. Measurement occurs only when the READ indicator is lit. Accuracy of the reading is $\pm 1\text{mm}$. Usable range is from 5mm to 35 mm.

- 4 Older animals with tough dead skin will require extra effort to obtain a good clean probe contact with the underlying live tissue. If necessary, remove the hair at the site. Then wet the skin with hot water. Next, apply light oil and allow it to penetrate the skin a minute or two before measuring. Young hogs, with softer skin, generally do not require much site preparation.

2. Why only measure the last rib?

Only site B (the last rib) will give an accurate reading of the three layers of fat. The lean indicator measures the skin and the first two or three layers of fat. Many hogs have the thin layer of third layer fat at an early age. As they increase in weight and age, the third layers become thicker and more wide spread over the body. Site B is the location normally used in metal ruler probing, and shows the best correlation with overall carcass yield. Because of the trapezius muscle at the shoulder (false lean), measurement at site C is more difficult and inconsistent. (True for both ruler and ultrasound devices.) Site A will give accurate measurement for two layers of fat but may not give an accurate measurement for three layers of fat due to the muscle tissue below the fat layers.

Values obtained using ultrasonic devices on live animals will generally measure 10% less than those measured at slaughter because measurement sites are somewhat different and the slaughtered carcass is disorted.

For comparison of backfat values between animals of different weights, the readings should be adjusted to a "standard" 100kg (220 lb.) hog. Tables of adjustment factors are readily available from most extension services. A convenient formula for the adjustment which is valid for a 100kg $\pm 15\text{kg}$ hog: **FAT READING x FACTOR**

$$\text{FACTOR} = 1.275 + ((.0033(\text{weight})) - (.0000605(\text{weight}^2)))$$

Probing should be done on gilts that are contemporaries and that are at or near market weight.

The lean-meater is powered by a NICAD rechargeable battery assembly which will typically last for 1,000 charges. Battery assemblies are specially designed for high reliability operation. For optimum performance replacements should be obtained from the factory. Battery is permanently connected to assure reliable operation under adverse conditions. Replacement should only be done by a qualified electronic service technician.

DANGER! Do not use the lean-meater when connected to charger.

WARNING! High voltage is present in unit. There are no user serviceable parts inside.

CAUTION! Do not use the lean-meater on human beings.

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