MS Multiscan

NL klantenservice
BE Klantenservice/Service Clientèle
DE Kundenservice
DK Kundeservice
IT Servizio clienti
ES Servicio al cliente
UK Customer Service
CA Customer Service
FR Service Clientèle
EX Customer Service Export

+31 (0)497 339 787
+32 (0)14-820713
+49 (0)2833-923630
+45 89884187
035-4490369
+34 931816433
+44 (0)01733592049
+1 866 995-7771
+33 (0)2 99 61 40 40
+31 (0)497 700 278

NL	MS Multiscan	
	Gebruiksaanwijzing.....	3
EN	MS Multiscan	
	Manual.....	25
DE	MS Multiscan	
	Gebrauchsanweisung	47
ES	MS Multiscan	
	Instrucciones para el uso	68
FR	MS Multiscan	
	Mode d'emploi.....	92

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EN MS Multiscan

1. Overview.....	26
1.1 Features.....	26
1.2 Setting up screen	26
1.3 Manufacturing Standards	26
2. System introduction.....	27
2.1 Outline	27
2.2 Schematic Diagram	28
2.3 Technical specification.....	29
2.4 Basic Principles	29
3. Installation.....	29
3.1 Environmental Requirements.....	29
3.2 Unpacking	30
3.3 Installation	30
3.4 Charging the Battery Pack.....	31
4. Operating panel + mouse.....	32
4.2 Operating Panel.....	32
4.4 Optional Mouse.....	34
5. Operating Instructions.....	34
5.1 Power On	34
5.1.1 Adjust Image Parameters.....	34
5.1.2 Frequency Setting	34
5.1.3 Gain Setting	35
5.1.4 Modify Brightness and Contrast	35
5.2 Settings note key.....	35
5.3 Image Process	36
5.3.1 Image Storage	37
5.3.2 Image Load	37
5.3.3 Image processing.....	37
5.4 Colour settings	38
5.5 Measurement of circumference, area and volume	38
5.5.1 Freehand method:	38
5.5.2 Ellipse Method:.....	39
5.6 Statistics	40
5.7 Distance measurement	40
5.8 Volume Measurement.....	41
5.8.1. 3-set volume measuring.....	41
5.8.2. 2-set volume measurement.....	42
5.9 Heart rate measurement (Only possible in "B/M" and "M" Mode	42
5.9.1 Heart Rate Measurement Distance method.....	42
5.9.2 Heart Rate Measurement Circumference and area method.....	43
5.10 Data input	43
5.11 photo shots + image printing	43
5.12 Power Off.....	43
6. Routine maintenance.....	43
6.1 Power Supply	43
6.2 Display Unit.....	43
6.3 Probe	44
6.4 Battery.....	44
7. Transportation + Storage	45
7.1 Environmental Requirements for Transportation + Storage.....	45
7.2 Transportation.....	45
7.3 Storage	45
8. Trouble shooting	45
8.1 Examination	45
8.2 Malfunction Solution	46

1. Overview

1.1 Features

The MS Multiscan uses the latest micro computer technology: digital scanning conversion (DSC), variable aperture, high dynamics, low noise, wide band receiver, log compression, dynamic filtering and border imaging. These features produce a clear, stable image with high resolution.

- Four display modes: B, B+B, B+M, M; Can realize image real time display, frozen, zoom; Hospital, Name, Age, Sex annotation; Measure Distance, circumference, area, volume, heart rate, gestation age Image gray scale 256 levels.
- Combined power supply mode of AC adapter and built-in Li-ion chargeable battery, 3 battery charging modes and the specialized brownout mode enables more lasting battery operation.
- 5 inches TFT-LCD display and programmable device (FPGA) and surface mounted technology (SMT) make this device compact and light in weight.
- PAL-D video output.
- With mouse mouthpiece.
- Jet moulding enclosure with hand-held structure makes it convenient for out diagnoses.
- The device consists of mainframe, probe and adapter.
- Standard configuration is S2/3.5MHz mechanical sector scanning probe, with S2/5.0MHz mechanical sector scanning probe, S2/3.5MHz and S2/5.0MHz waterproof mechanical sector scanning probes for option.
- The device is proved safe and effective via clinical validation.

1.2 Setting up screen

- Put the scanner on and press a random button (except reset).
- Press button 4 and then a white bar will light up on the right screen at the Mhz. This MUST be set to 5 (can be done with the use of button 3)
- Next setting takes place when you press button 4 again.
- Then the white bar will go to gain. This can also be set up with button 3 or 5
- Press button 4 again and make the next settin. (with button 3 or 5)

Advice is to always go the highest view (with button 3) and then 2 steps back (button 5) to get the ideal settings and image. Only the Mhz settings must always be on 5 Mhz, because the probe is a 5Mhz probe.

1.3 Manufacturing Standards

The MS Multiscan has been designed and manufactured in accordance with normal standards for electro-medical equipment. It carries the CE Mark.

The MS Multiscan displays the following symbols on the outside:

-  Type B equipment symbol
-  Class II equipment symbol
-  Caution! Please refer to the user manual
-  Mains power On symbol
-  Mains power Off symbol

2. System introduction

2.1 Outline


Figure 2.1 front view


Figure 2-2. Side View


Figure 2-3. Bottom Surface



Figure 2-4. Back Cover

2.2 Schematic Diagram

Figure 2-3. Schematic Circuit Diagram

2.3 Technical specification

Model		1609875 MS Multiscan digital	
Probes		Probe in standard	Optional probes
		S2/3.5MHzMechanical	S2/5.0MHzMechanical
Displayed Depth (mm)		192	192
Maximal detect depth		≥140	≥80
Resolution (mm)	Lateral	≤4 (depth≤80) ≤5 (80 < depth≤130)	≤3 (depth≤60)
	Axial	≤2 (depth≤80)	≤1 (depth≤60)
Blind area (mm)		≤8	≤8
Geometric position Precision %)	Horizontal	≤20	≤15
	Vertical	≤10	≤10
Monitor size (inch)		5.0 TFT-LCD	
Display modes		B、B+B、B+M、M	
Image gray scale		256 levels	
Depth (mm)		120 ~ 190	
Measure functions		Distance, circumference, area, volume, heart rate, gestation age	
Character display		Company name, Age, Sex	
Battery capacity		2200mAh	
Mainframe power		13W at non-charging operation / 25W at charging	
Adapter power		45W	
Mainframe		800g	

2.4 Basic Principles

The principles of B mode scanning are as follows: human/animal organs have differing densities and ultrasound speeds, i.e. their level of acoustic impedance varies (= average density x acoustic speed). The piezoelectric crystal (transducer) is triggered by an electrical pulse to produce an ultrasound wave of definite frequency. The ultrasound wave is transmitted into the body. An echo wave is formed when the ultrasound wave passes through the surface of different organs. These echoes, with different amplitudes, are received back via the same transmitting transducer. The receiving echoes (acoustic energy) are converted into electrical pulses. After the pulses are amplified, demodulated and DSC, etc., a standard video signal is formed and a sectional image of the organ is displayed on the monitor screen.

3 .Installation

3.1 Environmental Requirements

Temperature	: +10°C~+40°C
Humidity	: <80%
Atmospheric pressure	: 70kpa~106kpa
Power supply	: AC 230V±22V, 50Hz/60Hz or AC 110V±11V, 60Hz;
Display unit	: DC 14 V±0.5V

Keep the MS Muliscan away from powerful electrical fields, magnetic fields and high voltage equipment. Keep the screen away from direct sunlight. Shading the screen will improve image observation. Store the MS Muliscan in a well ventilated, dry and dustproof environment.

3.2 Unpacking

Carefully unpack the MS Muliscan and check the contents in accordance with the packing list, ensuring that there are no defective parts.

3.3 Installation

Check the power supply of the AC-adaptor and make sure it is in the proper range. Connect the power supply to the MS Muliscan. Making sure the power is switched off (power light off), connect the 5-pin probe connector to the MS Muliscan (the socket is located on the right hand side of the scanner). Switch on the MS Muliscan.

Caution

Do not connect or disconnect the plug for the power supply whilst the MS Muliscan is switched on. This could damage both the probe and the scanner. Avoid repeated unplugging of the probe once it has been connected to the scanner. This helps to prevent connection problems and damage to the socket and plug.


Figure 3-1. Attaching/Detaching Probe to Holder

- Attaching/detaching the probe to the holder:
The MS Muliscan has a holder to safely secure the probe. To attach, press the probe against the holder until it snaps into place. To detach, push the probe upwards and take it out (See figure 3-1)
- Install/Detaching the battery pack from the display unit:
Remove the battery cover by sliding the release catch, and then take out the battery pack (See figure 3-2)


Figure 3-2. install/Detaching the Battery Pack

3.4 Charging the Battery Pack

There are two ways to charge the battery pack

(1) Charging via Display Unit

Make sure that the battery pack is in the display unit. Connect the power cord of the "Input" port in the AC-adaptor to the power outlet, and then connect the power cord of the "Output DC14V/3A" port in the AC-adaptor to the "DC 14V" jack in the display unit. The "Charge" indicator (right beside the powerlight) on the front panel will light up and charging will begin.

When charging is complete, indicator light will go out.


(2) Charging via AC-Adaptor (optional)

Detach the battery pack from the display unit. Connect one end of the power cord to the "Charge DC 12.6V/1A" port on the AC-adaptor, and the other end to the "golden connector" on the battery pack. The "Charge" indicator on the AC-adaptor will light up and charging will begin.

Caution

Do not use the "Output DC14V/3A" port and the "Charge DC12.6V/1A" port on the AC-adaptor at the same time.

4. Operating panel + mouse

Picture 4-1. Screen Display

4.2 Operating Panel

Figure 4.2 shows the operating panel from the MS Multiscan


Figure 4-2. Operating Panel

FREEZE

Press it to shift between the states of freeze and real time.

Notes: **FREEZE** stands in the right lower corner, the image is frozen.

M/0 ~ **BM/9** are multifunction keys

4.3 Key functions

- During menu operation, they are used for selecting the sub-menu.
- While inputting Age and Time in the menu, they are used as numbers.(Further details are available)
- Use **M/0**, **ZOOM/1**, **+/2**, **B/7**, **BB/8**, **BM/9** to put names of the company and further details into the scanner .

As functions keys :

M/0 M Mode Scanning

In B/M mode, no matter in the state of real time or freeze, press it to enter B mode scanning, with an M displaying on the screen.

In non-B/M mode, press it to activate the menu of image processing and press **CLS** to exit.

ZOOM/1

In the state of real-time scanning, press it to alter the image magnification ratio. The depth is selectable between 120 and 190 in 8 levels.

+/2 Measure Mark

Press it in Freeze state, then “+ mark” will be displayed on the screen. Use direction keys or optional mouse to move it. (Refer to next chapter for further details)

▲3 ~ ▼6 Direction Keys

They are used for moving the cursor. **▲3**, **▼6** are used for altering the scanning depth in the real-time state, in “B” mode the same function as **ZOOM/1**. Press **▲3** to increase the depth, press **▼6** to decrease the depth, the current selected depth is displayed in the lower right-hand corner. While inputting names of company and further details, use **▲3**, **▼6** to turn the page up and down.

5<, **▶4** are used for activating the parameters on the right side, with activated parameters lighten. At this time, use **▲3**, **▼6** to alter.

B/7 B Mode Display


In the state of freeze or real time, press it to enter B mode (The default is single B).

BB/8 Double B Mode Display

In the state of freeze or real time, press it to enter double B mode. There are two B picture on the screen. One is “Freeze” picture, the other is “Real time” one; Repeat pressing, the picture can shift between “Freeze” and “Real time”. Press “Freeze”, then both of them show freeze picture.

BM/9 B and M Mode Display

In the state of freeze or real time, press it twice to enter B/M mode. Both pictures are displayed on the screen. The left is B real-time picture, and the right is M one. Repeat pressing, you can shift between B and M. On B picture, there is a vertical line composed by equidistant dots which is named sampling line.

Single B

BB

B M

M

Picture 4. 3.Four kinds of mode display

O.B. Obstetrics Table

Press it in freeze mode of “B” or “BB” to display the obstetrics menu. Press the number keys and further details are available in the following instructions. To exit, please press **CLS**.

DIST Measure Reference

Along with $\boxed{+/-}$, $\boxed{M/0}$ and direction keys, the measure of distance, circumference, area and histogram can be done. Further details are available in the following chapter.

NOTE

Press it, the note menu is displayed. Then press the number keys to enter the sub-menu. Refer to the following chapter for further details. Press $\boxed{CLS}$ to exit.

CLS Clean Screen

Press it to clean marks, notes and results. In menu status, press it to quit the menu.

RESET

Press it when the device crashes by accidents or operation mistakes, then the device can return to normal.

4.4 Optional Mouse

As a supplement to the keyboard, a mouse can make the measure quicker and more convenient. For this purpose always use a three button mouse.

The left, middle and right buttons of a three button mouse, all have their specific functions.

Left Button:

- The cursor is displayed when you press it and enter distance measure.

Middle Button:

- Press it in real-time or freeze mode, the function is the same as $\boxed{CLS}$.

Right Button:

- It can identify the start and end place of the cursor in the distance measure mode and shift between the two.
- It can identify the start position of the cursor in circumference and area measurement when using the free-hand method.
- The function is the same as $\boxed{DIST}$ in circumference and area measure by ellipse method.

5. Operating Instructions

5.1 Power On

Turn the MS Multiscan on by pressing the power switch on the right side of the display unit and the "power" indicator will light up. Pressing any key (except $\boxed{RESET}$) will enter scanning mode, after the start-up logo has appeared. Adjust Brightness, contrast and gain controls on the side of the mainframe to achieve the best visual display. Picture 5-1. Start-up Logo

Adjust brightness, contrast and gain controls on the side of the mainframe to achieve the best visual display.

5.1.1 Adjust Image Parameters

The parameters include frequency of the probe, gain, brightness and contrast. Press $\boxed{5}$ or $\boxed{4}$ in real-time state, one of them is lightened. Use $\boxed{\blacktriangle}$ and $\boxed{\blacktriangledown}$ to set the parameters and they are displayed in the upper right-hand corner.

5.1.2 Frequency Setting

Press $\boxed{5}$ or $\boxed{4}$ in real-time state to lighten frequency in the upper right-hand corner and use $\boxed{\blacktriangle}$ and $\boxed{\blacktriangledown}$ to adjust, with the range of 2.5MHz. 3.5MHz. 5.0MHz.

5.1.3 Gain Setting

Press 5 or 4 in real-time state to lighten gain in the upper right-hand corner and use 3 and 6 to adjust, in a range between 0 and 60.

5.1.4 Modify Brightness and Contrast

Press 5 or 4 to lighten “. ”, and use 3 and 6 to modify the brightness and contrast.

5.2 Settings note key

In freeze state, press NOTE, and the menu as shown in fig. 5.2.1. You can complete these functions:

V1.20

0.NAME
1.AGE
2.SEX
3.COMMENT
4.TIME
5.HOSP
6.LANGUAGE
7.ERASE
8.DEFAULT SET
9.BEEP-ON

V1.20 Software Version Number

- Press M/0, select “0.NAME” to input a name, as follow:

+

PLEASE ENTER NAME:

0-A 1-B 2-C 7-D 8-E 9-F

There are 26 letter keys and 1 space key available, with 3 and 6 you can turn the page up and down the select the necessary letter. Characters can be selected by pressing the correspondant number key, maximum amount of characters that can be inserted is 15 .

Press CLS to clear the wrong character. Press NOTE to confirm and exit after the text is inserted, or you can directly press NOTE to exit the input mode.

- Press ZOOM/1, select “1.AGE” to input the age or days of pregnancy, the maximum input number is 3 bits, the input screen is shown below:

PLEASE ENTER AGE:

Press NOTE to confirm and exit after inputting, or you can directly press NOTE to exit the input mode.

- Press +/2, select “2.SEX” to input sex of the animal, the input screen is shown below:

PLEASE ENTER SEX:
1.MALE 2.FEMALE

Use number key ZOOM/1 to select “1.MALE”, press +/2 to select “2.FEMALE”.

- Press 3 and select “3.COMMENT” to enter a image note.

White cursor will appear on the screen, together with 26 letter keys and space key at the bottom keys O.B. and DIST are used to turn page up and down in order to select the correspondent characters can be displayed when pressing the number key, with the maximum of 15 letters.

Press **CLS** to clean the wrong character. Press **NOTE** to confirm and exit after inputting, or you can directly press **NOTE** to exit from the input mode.

- Press **4**, and select "4.TIME" to modify time and date. The interface is as follow:

YY-MM-DD
HH-MM-SS

For example: 2009-3-22 9:35:30 are input as shown in the example below:

YY-MM-DD
090322
HH-MM-SS
093530

Press **NOTE** to confirm and exit after inputting, or you can directly press **NOTE** to exit from inputting.

- Press **5**, and select "5.HOSP" to input company name, as follow

PLEASE ENTER HOSP:

0-A 1-B 2-C 7-D 8-E 9-F

You can refer to name input for details how to switch and insert the company name correctly. There is a maximum of 18 characters in the name.

- Press **6**, and select "6.LANGUAGE". The interface is as

PLEASE ENTER LANGUAGE:
1.CHINESE 2.ENGLISH

Press **ZOOM/1** to select Chinese or **+2** to select English.

- Press **B/7**, and select "7.ERASE" to clear all image storage

ERASE ALL STORAGE?
1.YES 2.NO

Press **ZOOM/1** to confirm. During the process, "ERASING..." will appear in the upper left-hand corner and at the same time, other operations are unavailable. When it disappears, the storage is erased. Press **+2** to abandon and exit.

Attention:

Before the erasing is completed,(when "ERASING..." exists) , other operations are not allowed, they could damage the device.

- Press **BB/8** and select "8.DEFAULT SET" to restore the factory settings.
- Press **BM/9** and select "9.BEEP-ON" to switch on or switch off the beeping.

Note : When 0.NAME. 3.COMMENT. 5.HOSP are input , no number input is available. In note menu, press **NOTE** to exit directly.

5.3 Image Process

Press **M/0** in freeze state, the following menu enters the display . There are functions to save, load, image processing, change colour of the background and letters and to measure an area, volume and histogram.

V1.20: Software version number

V1.20
0.SAVE
1.SVLOAD
2.IMAGEPROC
3.COLOR
4.AREA-VOLM
5.HISTOGRAM

5.3.1 Image Storage

128 images can be stored permanently, also when the power is turned off for a longer period of time.

- Press **M/0** to display image process menu after a satisfactoring image is frozen, to store the image press **M/0**, the stored images automatically receive a number this number displayed in the upper left corner when the image is stored, example "SAVING.....05". When the storage is completed, the number disappears. Press **FREEZE** to return to real-time state.

When the memory is full(reaches 128 stored items), the following message will appear :

STORAGE IS FULL. ERASE NO.01 ?
1.YES 2.NO

If is the message states to erase No. 01 ; Press **ZOOM/1** to replace ; Press **+ / 2** to abandon the storage of the current image.

Notes :

- If the Memory for storage is full, and you want to store another item. A message will appear that an item should be erased in order to store the new image.

5.3.2 Image Load

Press **M/0** in the mode of real time or freeze to display image process menu. Then press **ZOOM/1** , the following will appear:

PLEASE ENTER STORAGE NO.:

Input the number of stored image according to the message, for example 01. Press **NOTE** after input, then No. 01 image is picked out. (If it is the wrong input, press **CLS** to clean the characters one by one and re-input) .01/128 is displayed in the lower left-hand corner on the screen. 01 is the number of current image and 128 stands for the storage capacity. At this time, press **▲ 3** or **▼ 6** to pick out images in other storage areas.


Press **FREEZE** to return to real-time state. To pick out other images, repeat the above procedure.

5.3.3 Image processing

Press **M/0** in real-time state to display the menu of image process, and then press **+ / 2** to select image processing. Repeat this to modify the smoothness and the image is displayed in the upper right-hand corner of the screen in real-time state. They are respectively: IM0. IM1. IM2. IM3.

Normal

Smoothed

5.4 Colour settings

Press **M/0** in real-time state to display the menu of image process, and press **▲3** to lighten colour on the right side. Use **◀5** or **▶4** to select modify items and press **▲3** and **▼6** to adjust the parameters.

IC0. IC1. IC2. IC3 are image coded colour, representing respectively grey, red, yellow and blue.

CC0. CC1 are colour of characters, representing grey and yellow.

BC0. BC1 are background colour, representing grey and blue.

5.5 Measurement of circumference, area and volume

Two methods are available.

- Press **M/0** in the state of freeze, the menu is displayed on the screen.
- Press **▶4** to select **4.AREA-VOLM**, and it is displayed as follow:

PLEASE ENTER:

1.FREEHAND 2.ELLIPSE

The first one is freehand method and the second is Ellipse method.

5.5.1 Freehand method:

Keyboard Operation

1. Press **ZOOM/1** to select Freehand method. The measure cursor appears on the screen. Use direction keys to move the cursor to the start of the examined spot.
2. Press **DIST**, and use direction keys to move the cursor along the fringe of the examined area to the end.
3. Press **DIST** again to finish the measuring.

If measuring should be continued, press **M/0** and **▶4** or directly press **+2**, and repeat Step 2-3. You can get 2 sets of values at most. The results are on the right side of the screen.

Picture 5.2. Illustration of Measuring circumference and Area(Free hand Method)

C1 and A1 are respectively circumference and area of the first set ;
C2 and A2 are respectively circumference and area of the second set ;
C1/C2 is the ratio of two circumferences ;
A1/A2 is the ratio of two areas.

Notes: There are limitations of measuring circumference and area by keyboard. For more conveniences, you can use the following Mouse Operation

Mouse Operation

1. Press **ZOOM/1** to select Freehand method. The cursor appears on the screen, use the mouse to move the cursor to the starting spot of the examined area

2. Press right button to move the cursor along the area that needs to be measured.
 3. Press right button again to complete the measuring of circumference.
- If measuring should be continued, press **M/0** and **▶4** or directly press **+2**, and repeat Step 2-3. You can get 2 sets of values at most. The results are mentioned on the right side of the screen. After the measuring, press Middle button to clear the screen.

5.5.2 Ellipse Method:

Keyboard Operation

Press **+2** to select Ellipse method. At this time an ellipse area appears which is called the examined area. Use direction keys to move this area. **DIST** is applied to shift among three functions of direction keys to adjust the size and angle.

1. Use direction keys to move the examined area along the displayed image.
2. Press **DIST**, and then direction keys to alter the size of the examined area. Press **▲3** and **▼6** to decrease or increase the area vertically, and then press **◀5** or **▶4** to decrease or increase the area horizontally;
3. Press **DIST** again, then use **◀5** or **▶4** to adjust the angle. press **◀5** to revolve anticlockwise the examined area and press **▶4** to revolve clockwise;
4. Press **DIST** again, the function of direction keys is to move the examined Area;
5. After the location, size, angle of the examined area are confirmed, the measuring can be done. If measuring should be continued, press **M/0** and **▶4** or directly press **+2**, and repeat Step 1 to 5. You can get 2 sets of values at most. The results are on the right side of the screen.

Picture 5.3. Illustration of Measuring circumference and Area (Ellipse Method)

C1 and A1 are respectively circumference and area of the first set ;
C2 and A2 are respectively circumference and area of the second set ;
C1/C2 is the ratio of two circumferences ;
A1/A2 is the ratio of two areas.

Mouse Operation

Press **+2** to select Ellipse Method. At this time an ellipse area appears which is called the examined area. Use direction keys to move this area. The right button is applied to shift among three functions of the mouse to adjust the size and angle;

1. Use the mouse to move the examined area to image display area;
2. Press the right button, and then move the mouse to alter the size of the examined area. Move the mouse left and right to decrease or increase the area vertically, and then move the mouse up and down to decrease or increase the area horizontally;
3. Press the right button again, then move the mouse left and right to revolve the examined area anticlockwise and revolve clockwise;

4. Press the right button again, the function of the mouse will be shifted into Move the Examined Area:
5. After the location, size, angle of the examined area are confirmed, the measuring can be done. If measuring should be continued, press **M/0** and **▶4** or directly press **+ / 2**, and repeat Step 1-5. You can get 2 sets of values at most. The results are shown on the right side of the screen. After measuring, press the middle button to clear the screen. Volume measuring is in the later chapter.

5.6 Statistics

- Press **M/0** in freeze mode, the menu is displayed on the screen.
- Press **5** to select **5.HISTOGRAM**, sampling window is displayed. Use direction keys or mouse to move to certain area you want to measure, and press **DIST** or the right button of the mouse to complete counting, the results are shown on the right side of the screen, as is shown in picture 5.4.

X axis stands for grey scale, and **y axis** stands for number.

PT stands for the total number of pixels in the rectangular window.

Gm stands for the grey scale of the curve at the peak of the y axis.

Pm stands for the number of pixels in Gm grey scale.

From the above figure, in the rectangular area, the total number of pixel dots is 10000. At gray scale 0, there are 962 dots, the most image pixel dots.

Picture 5.4. Illustration of the histogram statistics

- In the process, press **ZOOM/1** or **+ / 2** to reduce or enlarge the sampling window. Press **CLS** to exit.

5.7 Distance measurement

- Keyboard Operation :

1. Press **+ / 2** in freeze mode, the cursor is displayed on the screen.
 2. Use direction keys to move the cursor to the starting spot.
 3. Press **DIST** to confirm the starting spot for distance measuring.
 4. Press direction keys, another (second) cursor appears, move the cursor to the ending spot. Then press **ZOOM/1** to complete the measuring. (Note: Repeat pressing **DIST** to shift between the cursor of the starting spot and the ending spot)
- If distance measuring needs to be continued, you can repeat step 1-4, with the maximum of 4 sets of values. The results are displayed on the right side of the screen.

Picture 5.5. Illustration of Distance Measuring

The four sets are respectively D1, D2, D3, D4, in which
 $D1/D2$ is the ratio of D1, D2 ;
 $D3/D4$ is the ratio of D3, D4.

- Mouse Operation:

1. Press the left button and the cursor is displayed.
 2. Use the mouse to move the cursor to the starting spot.
 3. Press the right button for the starting spot.
 4. Use the mouse to move, and another cursor appears (cursor of the end point. Move it to the ending spot and press **ZOOM/1** to complete the measuring. (Note: Repeat pressing the right button can shift between the beginning spot and the end spot of the cursor)
- If distance measuring needs to be continued, you can repeat step 1-4, with the maximum of 4 sets of values. The results are displayed on the right side of the screen, as in Picture 4-5. Press the middle button to clean the screen after the measuring.

5.8 Volume Measurement

Two methods are available to measure the volume of a stored picture.

5.8.1. 3-set volume measuring

3 sets of distance values are measured by a 3 axis method, the result is obtained by calculation. The distance should be measured for three times before the volume can be measured, press **M/0** to obtain the result.

If the values are less than three sets, there will be no value displayed when you press **M/0** ; If you input four sets of values and then press **M/0** , the value displayed is the calculating result of the first three sets (D1, D2, D3) .

Procedures : (see picture 5.6 (kidney) as example)

1. Catch the cross and longitudinal sections of the kidney respectively and freeze them.
2. Measure the long axis and short axis of the cross section by means of distance measuring.
3. Measure the diameter of the longitudinal section by means of distance measuring.
4. Press **M/0** to complete the measuring, with the value of volume in "Vm1" on the right side, as follow :

Picture 5.6. Illustration of volume measuring (3 Axis Method)

5.8.2. 2-set volume measurement

Measure two sets of circumference and area by Ellipse Method and the result is obtained by calculating.

Procedures : (see picture 5.7 (Kidney) as example)

1. Catch the cross and longitudinal sections of the kidney and freeze.
2. Measure the circumference and area of cross and longitudinal sections.
3. The system will automatically complete the measuring, with the value of volume in "Vm1" on the right side, as follow :


Picture 5.7. Illustration of Volume Measuring
(Ellipse Method)

5.9 Heart rate measurement (Only possible in "B/M" and "M" Mode)

1. In B/M mode, freeze a desired image.
2. Measure the distance between the wave crests of two periods by means of the distance measuring method.
- 3 sets of values will be displayed in the lower right corner. The marks are respectively : Time T (unit : ms) . Heart Rate HR (unit : /m. Ejection Fraction EF unit : mm/s)

B/M Mode

Single M Mode

Picture 5.7. Illustration of Heart Rate Measuring

5.9.1 Heart Rate Measurement Distance method

"Left" click mouse to activate the cursor "+". Place the "+" at the start position of the measurement, "Right" click mouse to lock the cursor in place. Place the second "+" at the end point of the measurement (repeatedly "Right" click mouse to switch between the start and ending points). "Left" click mouse to get distance value. Repeat the above described method to obtain additional measurements.

5.9.2 Heart Rate Measurement Circumference and area method

Click the "Left" mouse button to activate the cursor "+". Move the "+" to the beginning of the area that needs to be measured, then Click the "middle" mouse button to lock the cursor in place at the starting point. Move the "+" to circle over the area that needs to be measured. "Middle" click mouse to get the value. Repeat the process to obtain additional measurements.

Note: In circumference and area measurement, when the trace is not totally closed, the starting point and the ending points are automatically connected to form a closed measurement trace.

5.10 Data input

To input data, "freeze frame" image, press **NOTE**. Shift cursor to "ID" and enter data – limit 6 numeral characters. After entering ID, the cursor will move to "AGE", enter data – limit 3 numeral characters. After entering age, the cursor will move to "SEX". Press **M/O** to input female symbol, or press **B M/9** to input male symbol.

Pressing **NOTE** again will erase all input data.

5.11 photo shots + image printing

- Photo Shots

Freeze frame image, then use Polaroid camera or digital camera to take photos.

- Image Printing

Connect "video out" of MS Multiscan to the video input on printer (connection cable sold Separately). Switch on power supply and operate printer in accordance with the manufacturer's instructions.

5.12 Power Off

Turn off the power supply switch. Disconnect the AC-adaptor if the scanner is not being used for a long time.

Note: Do not disconnect or re-connect the power supply until the scanner has shut down. Wait 2 or 3 minutes before restarting the scanner.

6. Routine maintenance

6.1 Power Supply

- Check the AC-adaptor power supply regularly. Do not use if voltage exceeds the recommended rating (AC220V±22V, 50Hz).
- Examine the AC-adaptor cable and probe cable regularly. Replace immediately if defective in any way.

6.2 Display Unit

- Operating conditions should be in accordance with the requirements detailed under "3.1 Environmental Requirements".
- Power off the display unit when cleaning. Use alcohol wipes to clean the case and an air duster to clean inner parts.

- Do not power on or power off too quickly - wait for 2 or 3 minutes before restarting the scanner.

If the MS Multiscan is not being used for a long period of time, re-pack it in its original case. Store the scanner in accordance with the conditions detailed under "7.1 Environmental Requirements for Transportation + Storage".

6.3 Probe

- The probe is expensive and can be easily damaged - avoid knocking or dropping the probe. After scanning, put the probe into its box and press **FREEZE** to keep it in the "frozen" mode.
- Only use "Ultrasound Transmission Gel" as the coupling medium for the probe.
- The probe itself is waterproof. It can be immersed in water in order to clean it after use. However, the connector must be kept clear of water during cleaning, and kept dry during use. Regularly check the body of the probe to ensure there are no cracks through which liquid might enter.
- Examine the probe frequently to ensure that it is completely filled with oil, without bubbles. Bubbles in the oil can occur, which will affect the image quality. When there are bubbles, please contact your supplying dealer in order to arrange for a repair to the probe.
- Cleaning, sterilization, and disinfection of probe: After scanning, use alcohol wipes to clean, sterilize and disinfect the probe, then return it to its box.
- Avoid unnecessary connection and disconnection of the probe to the display unit in order to minimize possible damage to plugs and sockets.

6.4 Battery

- The MS Multiscan is powered by a rechargeable polymer Li-ion battery.
- Please note, a new battery only reaches its optimum performance after 2 or 3 complete charge and discharge cycles.
- The battery can be charged and discharged hundreds of times but it will eventually wear out. When the operating time (scanning time and standby time) is noticeably shorter than normal, it is time to replace the battery.
- Use only batteries approved by the manufacturer of the MS Multiscan and only recharge batteries with an AC-adaptor approved by the manufacturer. Unplug the charger when not in use. Do not leave the battery connected to the charger for longer than 3 days - overcharging may shorten its lifetime. If left unused a fully charged battery will discharge itself over time.
- Temperature extremes can affect the ability of the battery to charge.
- Only use the battery for its intended purpose.
- Never use any charger or battery that is damaged.
- Do not short-circuit the battery. Accidental short circuiting can occur when a metallic object, ie: coin, clip or pen, comes into direct contact with either the + or - terminals (metal strips on the battery), for example when carrying a spare battery. Short circuiting the terminals may damage the battery or the connecting object.
- Leaving the battery in hot or cold places, such as in a closed car in summer or winter conditions, will reduce its capacity and lifetime. Always try to keep the battery between 15°C and 25°C (59°F and 77°F). Using an extremely hot or cold battery in the scanner may temporarily cause it to malfunction, even if the battery is fully charged. A battery's performance is particularly limited in temperatures well below freezing.
- Do not dispose of batteries in a fire!

- Dispose of batteries according to local recycling regulations. Do not dispose of batteries as household waste.
- Switch off the MS Multiscan before removing the battery.

7. Transportation + Storage

7.1 Environmental Requirements for Transportation + Storage

Environmental temperature	: - 5°C ~+55°C
Relative humidity	: <80%
Atmospheric pressure	: 80kpa~106kpa

7.2 Transportation

The MS Multiscan packaging is marked in accordance with the requirements of GB191 "Packing, Storage and Transportation Marks". The packing box contains shockproof protection for the MS Multiscan. The scanner packaging provides adequate protection for transportation by road, rail, sea or air. Exposure to rain, inversion and impact should be avoided.

The combination on the two locks of the carry case is set to "000" and locked in place with a plastic protector strip. To set the lock to a chosen combination, remove the plastic protector and slide the catch towards the middle of the case (the opposite way to opening it). With the slider held in this position, select the chosen combination by rotating the 3 wheels. Once the slider has been released, the chosen number will be the combination required to open the catch. Repeat the operation on each side of the case.

Ensure that these numbers are carefully noted.

7.3 Storage

- If the MS Multiscan has been stored for more than 6 months, remove the scanner from the carry case and recharge it for 4 hours. Return the scanner to the carry case in accordance with the instructions. Do not stack scanners on top of one another or store them directly against the floor, wall or ceiling.
- Store in a well ventilated, dry and dustproof environment. Avoid direct sunlight and caustic gas.

8. Trouble shooting

8.1 Examination

- Ensure the power supply is working normally. Make sure the power cable has been connected to the display unit correctly and inserted into the electric socket.
- Make sure the probe has been connected to the Multiscan correctly.

8.2 Malfunction Solution

No.	Malfunction symptom	Malfunction solution
1	The indicator light on the AC-adaptor is not on.	1. Check the power supply 2. Check the cable and plug of the AC-adaptor
2	The indicator light on the display unit is not on.	1. Check the cable and plug
3	Striations and snow appear on the screen.	1. Check the power supply of the AC-adaptor 2. Check the environment - make sure there is no interference from nearby electrical fields, magnetic fields or high-tension equipment 3. Make sure the power supply and probe plugs have been connected correctly
4	The image is not clear.	1. Adjust gain control 2. Adjust brightness and contrast controls

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