EN - Quickstart	2
ES - Quickstart	6
DE - Quickstart	10
FR - Quickstart	14
NL - Quickstart	18
DK - Quickstart	22
IT - Quickstart	26
EN - MS Smartscan Manual	30
PL - MS Smartscan Manual	68
EN - How to set the backfat reference Line on SmartScan	111

Wireless Vet Ultrasonic System SmartScan

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

1. Appearance


2. Connection of Host and tablet

2.1 Turn on tablet and read the manual of SmartScan


2.2 Press down and hold the Main Key for 7 seconds until the indicator turns blue to power on the machine.

The indicator will then flicker at a slow rate. Immediately the SmartScan will search for a matching WIFI signal pathway automatically.


2.3 Touch the icon of SmartScan APP on the screen to have the ultrasound interface on.

2.4 Successful Connection

---- When "Connected" appears on the bottom of the operating interface, it means that the smart device has been connected with SmartScan.

---- Press the Main Key to enter the real-time ultrasound scanning mode.

2.5 Unsuccessful Connection

---- If the following prompt pops up on the operating interface, it means that the smart device hasn't been connected with SmartScan. Reboot the SmartScan as instructed and click Retry until it's connected.

Hint

No ultrasound device found. Please turn on the device and try again

EXIT RETRY

---- If the same prompt appears on the operating interface after Retry several times, it indicates that the smart device is not installation-compatible and the SmartScan App is not properly installed.

3. Ultrasound Operation

Freeze/unfreeze

Tap this key to freeze/unfreeze the real-time image. Main Key of SmartScan can freeze/unfreeze the image as well.

Current Mode Selection

PT is for Pregnancy Test and BF is for Backfat measurement.

Setting

This App system offers user-defined pre-settings—Left/Right Hand and Auto Power-off Time.

Tap this key and the following window will pop up on the screen.

Load image

--- Tap this key and a dropdown menu will pop up on the screen. All saved images are listed in the menu;
--- Swipe the menu up and down to view and select the needed saved images.

frame_20170428_180807.jpg
frame_20170502_112812.jpg
frame_20170516_091122.jpg
frame_20170520_110824.jpg
frame_20170522_140832.jpg
frame_20170522_183805.jpg

The image is named as:

frame_Year/Month/Date_Hour/Minute/Second_.jpg

Save Image


To Freeze the image first before to tap this key to save it. The image is saved in jpg format on the tablet and its name is as:

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

Cine-loop bar

Once the images are frozen, the real-time images will be temporarily saved in the memory chip

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4 Ear-mark input. Ultrasound Operation

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
- 1) by manual
- 2) by reading ear-tag (RFID).
- 3) by scanning the barcode and QR

4.1 by manual

- Freeze the satisfied echo image first.
- Press the Save Key and the "Report" page will pop up on the screen:

- Click on the red blank ID space and a digit keyboard will appear at the bottom of the screen. Enter the corresponding ear mark and press "Next" (Blue Key) after confirmation.

- PT Result
Tap the optional box. Check the Pregnant box if the animal is pregnant. Leave the boxes blank if the animal is not pregnant.

Pregnant ☒ Not Pregnant ☐


- Remarks
--- Other information regarding to the animal or the image can be entered in the Remarks column. Enter words or sentences in English.
--- The Remarks column can be blank.
--- Press OK when all the information is entered.
- Check the report information and press the Ok key on the report page.
--- The report will then be saved in the smart device. Press the Cancel key to not save the report.
--- The path for saved Report: local RAM directory\MYJC\report\report_2017_05_23 09-23.pdf
--- The report files are named as:

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- View the images scanned in real-time and freeze the satisfactory one.
- Activate the Sub-Key and put the rear protecting part of the device close to the electronic ear tag worn by the animal.
---- Once the ear-tag is read, a beep from the SmartScan will be heard and the ear-tag number will be instantly displayed in the ID area on the left upper corner of the screen.
---- The ear-tag must be read within 10 seconds after Sub-Key is activated.
---- The ear-tag can be automatically read as shown below:
- PT Result
- Remark s
- Check the report information
(Complete Step 3 to 5 as instructed in by Manual)

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

If your machine is Version A and B, this function is available.


1) Press the PT/BF key to switch system to BF mode. The characters BF will be display on the right upper corner of the screen in steady of PT.

2) Place the probe on P2 point of the animal.


Obtain a correct ultrasound image and then freeze it.

3) After the image is frozen, to press on the upper area of the screen gently, a green line will be appeared automatically on the image. The end of the green dotted line is the lower end of the BF line. The length of the green line represents the thickness of BF.

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Strongly Suggestion:


In order to obtain the better image of the fat, you are strongly suggested to have the probe put on P2 at the right direction:

6. Pregnant Report

Whether in manual entry of ear marks, automatic reading of electronic ear tags, barcode or QR code, the users must check the test results carefully and enter the information truthfully. After the information is checked, tap the Ok key on the report page and a PT report of the tested animal will be automatically generated in the MYJC folder of the local root directory of the tablet.

Save path: RAM directory\MYJC\report\
Save format: .pdf

7. Database

In the RAM\MYJC\report\ folder, a database will be automatically generated as well. In the excel format, the file is named report.xls.

- 1) It's a truthful record of each test conclusion after the ear mark is entered;
- 2) The database is automatically generated and the contents can't be altered;

3) If the test was not conducted this time, the corresponding item will be blank.

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
Search path:

All saved files can be sent or transmitted through the following functions of the tablet:

- 1) Bluetooth
- 2) WIFI
- 3) USB (digital cable is packed in the white small carton box)

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MADE BY  solown

Contents

Foreword.....	1
About this Manual.....	1
Manufacturer's Statement.....	1
Manufacturer's Warranty.....	1
Notes.....	2
Classification.....	3
Product Labels.....	3
General Instructions for Device Operation.....	3
General Safety Information.....	4
Scope of Application.....	5
Contraindications.....	5
1 Product Information.....	6
1.1 Appearance and Dimensions:.....	6
1.2 Product Classification.....	7
1.3 SmartScan App.....	7
1.4 Related Parameters.....	7
SmartScan Ultrasound Technical Specifications.....	7
Battery Parameters.....	8
1.5 Hardware.....	8
1.6 Software.....	9
Network Safety.....	9
System Requirements.....	9
2 Functions and Operations.....	11
2.1 Function Keys.....	11
2.1.1 SmartScan Function Keys.....	11
2.1.2 Smart DeviceInterface.....	12
2.2 Connection with Smart Device.....	13
2.2.1 Get an Installation-Compatible Smart Device.....	13
2.2.2 Installation of SmartScan App.....	14
2.2.3 Enter Ultrasound Operating Interface.....	14
2.2.4 Power on the SmartScan.....	15
2.2.5 Successful Connection.....	15
2.2.6 Unsuccessful Connection.....	15
2.3 Ultrasound Diagnosis Operation.....	16
2.3.1 Freeze.....	16
2.3.2 Mode Selection.....	16
2.3.3 Gain and Depth Adjustment.....	17
2.3.4 Settings.....	17
2.3.5 Load Image.....	17
2.3.6 Save Image.....	18
2.3.7 Save Video.....	18

2.3.8 Cine-loop.....	18
2.3.9 Cine-loop Bar.....	19
3 Ear-mark Input.....	20
3.1 Input by manual.....	20
3.2 Reading Ear-tag (RFID).....	22
3.3 Barcode and QR Code Scan Entry.....	23
3.4 Report Generation.....	24
3.5 Database.....	24
3.5 Transmission of Saved Files.....	25
4 Backfat Measurement.....	26
4.1 Steps for backfat Measurement.....	26
4.2 Save the result of backfat measurement.....	28
5 Storage & Backup of the saved files.....	29
5.1 File Saving Formats.....	29
5.2 File Saving Path and Folder.....	29
5.3 Data Backup.....	29
5.4 The safest way to backup.....	30
6 Power Charging.....	31
Battery Replacement:.....	32
7 Maintenance	32
7.1 Routine Check and Maintenance.....	32
7.2 Cleaning and Disinfection.....	33

Foreword

About this Manual

This Manual is included with the purchase of SONOWIN's SmartScan. Any form of plagiarism, duplication or clipping of this Manual is prohibited without written consent from SONOWIN.

This Manual contains the following contents:

1. Product Information
2. Functions and Operations
3. Ear Mark/Ear Tag Entry
4. Back Fat Measurement
5. File Saving and Backup
6. Battery Charging
7. Maintenance Instructions

Product Model of this Manual: SmartScan

Version: SM. V1.00e

Manufacturer's Statement

- 1) This Manual and the photos and illustrations contained herein are protected by international copyright laws, and all rights are reserved by Shenzhen SONOWIN Science Technology Inc.
- 2) Without the prior written consent from the manufacturer, this document should not be photocopied, duplicated or translated into other languages.
- 3) The contents of this Manual may be revised by Shenzhen SONOWIN Science Technology Inc. without prior notice.
- 4) Certain photos contained herein are schematic diagrams for reference only. In case of inconsistency with the product, the product shall prevail.
- 5) SONOWIN is the registered trademark of Shenzhen Sonowin Science Technology Inc. Unauthorized use of the trademark will be legally prosecuted.

Manufacturer's Warranty

- 1) Shenzhen SONOWIN Science Technology Inc. will only be responsible for the safety, reliability and performance of the device, provided that the device is dismantled and serviced by the designated professionals of Shenzhen SONOWIN Science Technology Inc. and used in strict accordance with this User's Manual.
- 2) Shenzhen SONOWIN Science Technology Inc. undertakes that the warranty of the device is two (2) year starting from the date of purchase and delivery, and that the brand new device is free of material and process problems. Within the warranty period, Shenzhen SONOWIN Science Technology Inc. will repair the device and replace parts damaged not due to human use free of charge. However, devices with surface damage will not be repaired or replaced.

- 3) This warranty is only applicable to failures resulting from operation in accordance with the designated conditions within the specified range of application set forth herein.
- 4) Any loss or damage of the device caused by external factors like lightning strike, earthquake, theft, improper use or abuse and modification is not covered by this warranty.
- 5) Shenzhen SONOWIN Science Technology Inc. will not be held responsible for damage arising from other devices or unauthorized connection with other devices.
- 6) Shenzhen SONOWIN Science Technology Inc. will not be held responsible for the loss, damage or injury arising from delay of service requests.
- 7) In case of product problems during the warranty period, please contact Shenzhen SONOWIN Science Technology Inc., and state the device model, serial number, date of purchase and nature of problem.

Notes

To ensure operational safety and long-term stable performance of the device, please read this Manual carefully before use to fully understand the functions, operations and maintenance tips for the device. Attention must be paid to contents marked Warning, Caution, or Notes herein.

◆ Misoperation or noncompliance with the instructions or the manufacturer or the agents thereof may result in damage of the device or personal injury.

◆ The following terms used in this Manual are to stress particular information:

Warning — Ignoring this information may lead to severe personal injury, death or property loss.

- a. Do not sell or use this product in a strong nuclear magnetic field.
- b. The head of the device is an ultrasound probe. Don't hit, drop or exert excessive pressure on it. Always wear the wristband during use.
- c. Do not dismantle the product or the wiring of any part.
- d. Do not use the device in case of damage of the wires or charger. Return the damaged parts to SONOWIN for replacement.
- e. The device contains a built-in lithium ion battery. Do not use it in an inflammable environment to prevent explosion.

Cautions—Ignoring this information may lead to minor personal injury or property loss.

- a. To reduce the electromagnetic interference from surrounding devices, please use this device in conditions with the least possible instantaneous power supply, mechanical impact, vibration, thermal radiation, optical radiation and ion radiation from surrounding devices.
- b. Use anti-static spray on the carpet and felt, and use an anti-static mat to avoid damage from static electricity.
- c. Despite of the IPX7 waterproofing rating of this device, please store or place it out of direct sunlight and away from moisture.
- d. Please use a medical ultrasound coupling agent. Non-compatible coupling agents will not only affect the ultrasound image, but also cause skin discomfort of the user or animal.
- e. Do not expose the device to direct heat sources like strong sunlight or local heat sources to prevent aging or decrease in sensitivity.
- f. This product is a precision device and thus should be handled with care. Please contact SONOWIN for repair in case of any damage.
- g. Do not immerse the product for an extended period. If the product is immersed, please do not start it up and use it; contact SONOWIN for maintenance.
- h. The SmartScan App is the intangible asset of SONOWIN and should not be transferred or traded.

Notes — Important installation, operational or maintenance information, but not dangerous. Warnings against danger are not given in Notes.

Classification

Based on Electric Shock Resistance Type: Type I externally powered device

Based on Electric Shock Resistance Level: BF application part

Based on Seepage Protection: IPX7; please refer to relevant manuals for smart phones or tablet PCs.

Product Labels

The following labels are marked on the products of SONOWIN and the accessories, spare parts and even packages thereof, and have the meanings as described below respectively:

Lable	Name	Lable	Name
	Manufacturer		易碎标识
	Model Name		Waterproof
	Product Serial Number		Storage Temperature
	QC Pass Certificate		Stored Upright
	Warning		Certification Label
	EU Compliance		

General Instructions for Device Operation

During Operation

1. Please have the wrist **strap on your hands during** using this device to avoid any falling.
2. Keep an eye on the changes of Main Key indicator to make sure the device is under wifi connection.

3. Do not drop or hit the device during use.
4. Wait at least 2 to 3 minutes after shutdown before restarting the device.
5. In case of abnormality during ascan, immediately stop and shut down the device.
6. Do not put excessive upward or downward pressure on the device during use to prevent possible damage to the sound head oiler.
7. Always check the battery level and charge the battery in a timely manner when necessary.
8. The device can be used during charging, and extra attention should be paid to safety.

After Operation

1. Cut off the power.
2. Check whether the waterproof rubber plug of the charging base is properly closed before washing the device.
3. Gently wipe the probe oiler with a piece of wet soft cloth. Try to reduce friction damage on the sound head oiler.

General Safety Information

■ The following safety measures must be taken to ensure the safety of the operator and patient and the reliability of the device:

1. The device should be used by or under the direction of a qualified operator.
2. Do not change the parameters of the device. If change of parameters is required, please contact a local dealer or SONOWIN.
3. The device is tuned to best performance before it's delivered. Use it as instructed by the Manual. Do not adjust any preset control or switch.
4. In case of failure of the device, immediately shut down the device and contact a local dealer or SONOWIN.
5. To upload the product information to the internet, please contact a local dealer or SONOWIN for value-added service.
6. Operating, Storage and Transport Environment of the Device:
 - Avoid violent vibration under normal operation and keep the temperature, humidity and power supply within the following range:
 - 1) AmbientTemperature Range: -20°C-+60 °C
 - 2) Relative Humidity Range: 30%-80%
 - 3) Power Supply Range: a.c. 100V-240V, 50/60Hz
 - Avoid exposure to the rain or snow and mechanical collision during transport. Do not load or transport the device with corrosive substances.
 - Since the host contains a built-in battery, the device should be stored in a dry place:
 - 1) AmbientTemperature: -20oC-+50 oC (for less than 1 month of storage)
AmbientTemperature: -20oC-+45 oC (for 1 to3 months of storage)
AmbientTemperature: -20oC-+20oC (for less than 1 year of storage)
 - 2) Relative Humidity: Less than 80% (20oC)
 - 3) The indoor storage environment should be protectedfrom strong sunlight and other corrosive gasses, and well-ventilated.
7. To reduce the risks during use, please use the power cordthat comes with the device and standard grounded grid power supply.
8. SONOWIN will assume no liability for risks arising from unauthorized dismantling or modification of the device by the user.
9. Turn off the power if the device is not to be used for a long period of time.
10. Do not press the probe too hard on the patient during the diagnosis to avoidunnecessary injury.

Scope of Application

In vitro ultrasound diagnostic examination for animals.

Contraindications

The device is not applicable to examination of gas-containing organs and tissues of animals, e.g. lungs.


1 Product Information

This product is a general handheld software-controlled veterinary ultrasoundscanner. It obtains real-time ultrasound signals through WIFI and reflects them on Android devices like smart phones or tablet PCs.

The product features a full touch screen on the operating interface of the smart device.


This product is composed of a built-in battery, an ultrasound transducer, a multi-channel all digit synthesizer, WIFI components, RFID module and coil. A battery adapter and a lockable wristband come with the product. If a tablet PC is purchased with the product, a cover for the tablet PC and a neckband will be additionally provided.

1.1 Appearance and Dimensions:

Host (Probe)

User Interface/Screen on the Smart Device

Front

Back

Left Side View

Right Side View

Item	Length (mm)	Width (mm)	Height (mm)	NW
Host (Probe)	190	Max. 50	Max. 50	285g
Smart Tablet PC (box)	21	13	12	350g
Package (Host + Tablet PC+ Support Panel + Mini-holder+ Wristband+ Tablet PC Accessories+ Neckband)	240	180	140	2.7KG

1.2 Product Classification

High Version (A) — PT+ BF+RFID

Standard Version (B) — PT+ BF

Low Version (C) — PT only

13. SmartScan App

Software Version: MYJC-V1.00

1.4 Related Parameters

SmartScan Ultrasound Technical Specifications

Technology	Full digital beam former
Scan Method	Mechanical sector scan (with negative pressure controlling)
Probe Frequency	3.5MHz
Display Mode	B
Gray Scale	256-Level
Depth Control	100-180mm, adjustable
Gain Control	90-120dB, adjustable
Image Procession	Endge enhancement, frame rate related
Cine Loop	Max. 128 frames
Saved Image	.jpg format
Saved Video	.mp4 format
Pregnancy Report	.pdf format
Testing Data	Excel format for pregnancy and backfat Conclusion saved automatically

Backfat Thickness	Auto identification, auto measurement and auto saving the results
Ear Number Input	1) by manual
	2) by barcode scanning
	3) by QR scanning
	4) by RFID reading, FDX & HDX both available (optional for Version A only)
Power Consumption	0.65mA, 3.6V
Adapter Input	VAC100-240, 50/60Hz, 5V, 2A
Working Time of Battery	More than 3 hours (under continuous real-time scanning)
Package Size	240*180*140 mm, carton box
Host Net Weight	285g
Package Gross Weight	2.7Kg

Battery Parameters

Battery Manufacturer	Shenzhen Syren Energy Technology Co., Ltd.
Battery Type	SANYO Lithium Ion Battery
Battery Model	NCR18650B-H00HA
Battery Management:	JEITA compatible charger, packaged fuel gauge and protective circuit, battery balance and temperature monitoring
Charging Input:	VAC100-240, 50/60Hz, 5V, 2A
Battery Capacity:	3400mAh

1.5 Hardware

Purchase

- The life of the host (probe) is about 6 years.
- Please contact your local exclusive distributor MS Schippers (<http://www.schippers.eu>) to purchase other value-added services or optional accessories.
- Please contact the local distributor in your market.

Warranty

- The warranty period of respective parts of the device is as follows:

Item	Warranty Period	Remarks
Host (Probe)	2 Years	
Smart Android Tablet		Reference its warranty advice
Other Accessories	1 Year	

1.6 Software

- ◆ The icon of the SmartScan App is:

- ◆ If a smart device is presented during the purchase, the SmartScan App will be properly installed and ready to use after boot;
- ◆ If no smart device is presented during the purchase, the SmartScan App should be provided to the buyer or a local dealer for installation as required.
- ◆ The current version of the SmartScan App is MYJC-V1.0. Please contact [local exclusive distributor MS Schippers \(http://www.schippers.eu\)](http://www.schippers.eu) for software updating.

Network Safety

To ensure the stable and safe connection between a smart device and SmartScan, the SmartScan App adopts the self-configured WiFi connection mode for automatic encryption and connection.

In places with concentrated WiFi signals, the self-configured WiFi signals will be inevitably affected and the SmartScan may experience frame rate lowering. There is no need to worry about this normal phenomenon.

Suggestion:

Do not use the device in places with concentrated WiFi signals.

Once the WiFi is disconnected, the smart system will automatically shut down after a period of continuous monitoring to enhance the working hours of the battery.

System Requirements

The smart device on which the SmartScan App runs must meet the following minimum requirements:

Technical Specifications:

- Supports BLE (Bluetooth) v4.0+
- Supports WIFI 802.11n and WIFI direct use
- At least 16G hard disk
- At least 2G RAM

Operating System:

- Android System 6.5 (Only compatible with systems below 6.5)

Processor:

- Dual core CPU
- ARM-based CPU architecture (Android device)

Display:

- IPS, Resolution (Pixels): 1280*800
- Support OpenGL ES 2.0

Supportable Smart Devices:

- Android System 6.5 (Only compatible with systems below 6.5)
- Devices with WIFI 802.11n
- Devices with BLE (Bluetooth) 4.1


- The SmartScan App must be properly installed on the smart device.
- If the smart device has insufficient memory space, certain functions of the device can't be enabled.
- In consideration of the openness and complexity of Android System, SONOWIN will not provide customers with the App for self installation.

2 Functions and Operations

2.1 Function Keys

2.1.1 SmartScan Function Keys

The main key has the following three functions:

- Power on/off
- wifi enable
- Image frozen/defrozen


The sub-key is for switch on/off reading ear-tag.

- it is available on premium version A only.
- It works only when image is frozen.

Key Name	Function	Operations
 Main Key	Power On/Off	Press down and hold the key with the thumb for about 7 seconds to power on the SmartScan. The main key will turn in blue and its indicator will flicker at a slow rate. Immediately the SmartScan will automatically search for matching WIFI. During use, press and hold the key for 7 seconds to power off the SmartScan. The indicator of the MainKey will then be turned off and the power disconnected.
	wifi Enable	The WIFI will be enabled when the device is powered on. Once the WIFI is connected, the blue indicator of the PrimaryKey will flicker at a fast rate. "Connected" will flash on the screen. If the WIFI is not connected, the blue indicator of the PrimaryKey will continue to flicker at a slow rate. After flickering for 1 minute, the device will shut down automatically.
	Real-time/ Freeze	After the WIFI is connected, press the Real-time/Freeze Key to enter the real-time scan mode; Then press it to freeze/unfreeze the image. It works the same way as the Freeze Key on the smart interface.

 Sub-Key	RFID On/Off	This key is available for the premium version A only. Notes: It is available to be used when the image is frozen only. --- To freeze the image before pressing this key. Then, to have the rear part of the device close to the ear-tag to read out its number within 10 seconds. The reading distance is around 6-10 cm.
--	----------------	---

2.1.2 Smart Device Interface

Icons	Description
Freeze	Tap the icon Freeze during the real-time scan to freeze the image, tap it again to unfreeze the image. (The main Key of SmartScan also has the same function)
PT/BF	PT is short for Pregnancy Test . BF is for Backfat measurement. SmartScan comes three versions in all. Both the premium version A and standard version B have this key. The power-on default setting is PT mode which can be seen in the upper right corner. Press it to switch to BF mode which can be seen in the upper right corner.
Setting	Tap this icon and two settings will pop up: (1) Left/Right Hand Setting (The default setting is Left Hand); (2) Screensaver Time Setting (1-20 seconds available to choose from).
Load Image	Tap this icon to load the saved images.
Save Image	Tap this icon to save the frozen image on the smart device.

Save Video	Tap this icon to save the cine-loop images. Reminder: This function can only be enabled when the smart device has 2G available memory space.
Cine Loop	After freezing the image, tap this key to cine-loop the images . The maximum is 128 frames ---- Tap this key again to pause the playback. ---- Swipe the screen left or right to play back the image frame by frame; ---- The temporary saving number of the current image is prompted as 023/ 128.
Cine-loop Bar	After the images are frozen, the number onthe progress bar indicates the maximum frame number of the real-time scanned images before freezing. Roll the slider and the frame number of the currently displayed image will change accordingly.
Battery Power	This icon indicates the battery power of SmartScan; --- When the battery is fully charged, there will be 4 bars. --- Please to charge the SmartScan when there is only one bar left.
Current Mode	This icon displays the current mode. PT -- Pregnancy Test; BF -- Back Fat measurement Noticed: This icon does not work in Version C
Date/Time	Date: Year/Month/Date (Year is shown in 4 digits and month and date in 2 digits; “0” will be prefixed if there is only 1 digit.) Time: 24-hour base
Gain	This icon displays the current gain information. Swipe the right part of the screen up and down to adjust the current gain.
Depth	This icon displays the current depth information. Swipe the left part of the screen up and down to adjust the current depth.
ID	ID means the ear mark of the animal. ---- The ear mark could be as ear-tag, or barcode, or QR or any character input by manual. So there are four (4) difference way tobe entered: 1) by reading ear-tag with RFID tech ---- Only available in SmartScan A 2) by inputting manual ---- Available in the whole SmartScan series. 3) by scanning the barcode --- Available in the whole SmartScan series. 4) by scanning QR code --- Available in the whole SmartScan series.

2.2 Connection with Smart Device

2.2.1 Get an Installation-Compatible Smart Device

---- If a tablet PC is purchased with SmartScan together, it will have been installed with the SmartScan App;

---- If a tablet PC is not purchased with SmartScan, buyers are required to get an installation-compatible smart device and request the SmartScan App from SONOWIN or a local dealer for self-installation.


---- The interface of the ultrasound system can be switched into difference language when the smart bablet is switched to.

Suggestion:

- Buyers are advised to purchase the smart device together with SmartScan to avoid installation problems.
- The smart device offered is Huawei Honor 2, 8 inch screen.


2.2.2 Installation of SmartScan App

The SmartScan App can be easily installed. Click Install for automatic installation. The orange icon of SmartScan will be automatically displayed on the screen after installation.


2.2.3 Enter Ultrasound Operating Interface

Click the SmartScan icon to enter the Ultrasound Operating Interface automatically like below:


2.2.4 Power on the SmartScan

- Press and hold the Main Key for 7 seconds until the indicator turns blue to power on the SmartScan;
- The indicator will then flicker at a slow rate. Immediately the SmartScan will automatically search for a matching WIFI signal pathway.
- Check the bottom of the screen on the smart device.


2.2.5 Successful Connection

- When “Connected” appears on the bottom of the operating interface, it means that the smart device has been connected with SmartScan.
- Press the Main Key to enter the real-time ultrasound scan mode.

2.2.6 Unsuccessful Connection

- If the following prompt pops up on the operating interface, it means that the smart device hasn't been connected with SmartScan. Reboot the SmartScan as instructed and click Retry until it's connected.
- If the same prompt appears on the operating interface after Retry several times, it indicates that the smart device is not installation-compatible and the SmartScan App is not properly installed.

2.3 Ultrasound Diagnosis Operation

Press the PrimaryKey after the SmartScan is connected to WIFI of the smart device to immediately enter the real-time ultrasound scan mode.

2.3.1 Freeze

Tap this key to freeze/unfreeze the real-time image.
The Main Key of SmartScan can freeze/unfreeze the image as well.

2.3.2 Mode Selection

This key is only available in the premium version A and standard version B.

PT is for Pregnancy Tes, and BF is for Backfat measurement.

The power-on default mode is PT mode which can be seen in the upper right corner. Press the PT/BF key to switch to BF mode which can be seen in the upper right corner.

Preset Parameters on PT mode:

Depth: 160 mm
Gain: 100 dB


Preset Parameters on BF mode:

Depth: 100 mm
Gain: 90 dB

2.3.3 Gain and Depth Adjustment


Swipe the left and right parts of the screen up and down during the scan for gain and depth adjustment.

<u>Left: Depth</u>	<u>Right: Gain</u>
Adjustable Range: 100-180mm	Adjustable Range: 100~120dB
Increase/Decrease by 10mm each swipe	Increase/Decrease by 1dB each swipe

2.3.4 Settings

This App system offers user-defined pre-settings—Left/Right Hand and Auto Power-off Time. Tap this key and the following window will pop up on the screen.


Tap on “Left” to switch to “Right” for right hand-held smart device operation.

Choose from 1 to 12 for the Auto Power-off Time. The default value is 9 seconds.

Press Ok or Cancel to confirm the settings.

2.3.5 Load Image

Tap this key and a dropdown menu will pop up on the screen. All saved images are listed in the menu;

Swipe the menu up and down to view all saved images;

Images are named as:

frame_Year/Month/Date_Hour/Minute/Second_.jpg

Click the selected file and the image will be immediately displayed on the screen.

2.3.6 Save Image

To Freeze the image first before to tap this key to save it on the smart device.

The path for the saved images: local RAM directory\myjc\Image\frame_20170522_140988_.jpg

The format of the saved images is in .jpg format.

The saved images names as:

Frame_Year_Month/Date_Hour/Minute/Second_.jpg

2.3.7 Save Video

To Freeze the image before tap this key to automatically save it on the smart device.

Videosave path: local RAM directory\myjc\Image\record_20170505_153443_.mp4

The saving format of videos is mp4.

The video files are named as:

record_Year/Month/Date_Hour/Minute/Second_.mp4

Reminder:


Videos can only be saved when the smart device has 2G available memory space.

2.3.8 Cine-loop

After freezing images, tap this key to play back the images scanned in real-time before freezing one by one; 128 frames maximum for playback;


Tap this key again to pause the playback. Swipe the screen left or right to play back the image frame by frame;

The temporary save number of the current image is dynamically prompted as 023/128 behind the progress bar.

2.3.9 Cine-loop Bar

Once the images are frozen, the real-time images will be temporarily saved in the memory chip; The maximum temporary frames will be automatically displayed behind the cine-loop bar. Roll the slider for quick playback.

3 Ear-mark Input

It would not be the practically meaningful if the ear-mar (or ID) of the tested animal does not be input whenever the frozen images or backfat measuring conclusions is saved.

The saved image will saved on an indpenent “Report” page. The conclusion of diagnosie and the measuring results will be saved in an Excel datababes. Both the Report and the data can be transferred from the smart table to any other device and PC.

The ear-mark can be entered in the following four ways:

- 1) by mannual
- 2) by readinnng ear-tag automattically .
- 3) by scaning the barcode automatically
- 4) by scanning QR automatically

Important:

The images must be frozen before they can be entered.

SmartScan comes in the following difference versions:


- (1) high version A
- (2) standard version B
- (3) basic version C

Version A is with RFID for both FDX and HDX.


3.1 Input by manual

To save the scanned image and its diagnosis conclusion and send it to the management system, the ear mark of the tested animal must be entered as follows:

- 1) View the images scanned in real-time and freeze the satisfactory one.
- 2) Press the Save Key and the “Report” page will pop up on the screen:

3) Click on the red blank ID space and a digit keyboard will appear at the bottom of the screen. Enter the corresponding ear mark and press “Next” (Blue Key) after confirmation.

(The animal ID number can be a random combination of digits and characters with no limitation on characters.)

4) PT Result

Tap the optional box. Check the Pregnant box if the animal is pregnant. Leave the boxes blank if the animal is not pregnant.

Pregnant ☒ Not Pregnant ☐

5) Remarks

Other information regarding to the animal or the image can be entered in the Remarks column. Enter words or sentences in Chinese.

The Remarks column can be left blank.

Press OK when all the information is entered.

(Words or sentences can be entered in the Remarks column based on the actual conditions of the farm.)

- 6) Check the report information and press the Ok key on the report page.
The report will then be saved in the smart device. Press the Cancel key to not save the report.

The path for saved Report: local RAM directory\MYJC\report\report_2017_05_23 09-23.pdf
The report files are named as:

record_Year/Month/Date_Hour/Minute.pdf


- 7) See 3.5 Single Report for more report details.

3.2 Reading Ear-tag (RFID)

SmartScan in Verions A can read the ear-tag by switching on the sub-key RFID instead of manual entry.

Specifically:

- 1) View the images scanned in real-time and freeze the satisfactory one.
- 2) Activate the Sub-Key and put the rear protecting part of the device close to the electronic ear tag worn by the animal.
 - Once the ear-tag is read, a beep from the SmartScan will be heard and the ear-tag number will be instantly displayed in the ID area on the left upper corner of the screen.
 - The ear-tag must be read within 10 seconds after Sub-Key is activated.\
 - The ear-tag can be automatically read as shown below:
- 3) PT Result
- 4) Remark s
- 5) Check the report information
(Complete Step 3 to 5 as instructed in by Manual)


3.3 Barcode and QR Code Scan Entry

Barcodes in the following two formats can be scanned by SmartScan:

Barcode


QR Code:

LYHHXM150105184

How to Scan:

Tap the Camera icon and the tablet PC will automatically start the scan process; point the scanner of the smart device at the barcode and QR code to be scanned. A beep from the smart device within 1 to 3 seconds indicates a successful scan.

3.4 Report Generation

Whether in manual entry of ear marks, automatic reading of electronic ear tags, barcode or QR code, the users must check the test results carefully and enter the information truthfully.

After the information is checked, tap the Ok key on the report page and a PT report of the tested animal will be automatically generated in the MYJC folder of the local root directory of the smart device.

Save path: RAM directory\MYJC\report\

Save format: .pdf

Report

Date: May 23, 2017
ID: MYHC01A170603001
Pregnant: Yes

Remarks:
First pregnancy

3.5 Database

In the RAM\MYJC\report\ folder, a database will be automatically generated as well. In the excel format , the file is named report.xls.

- 1) It's a truthful record of each test conclusion after the ear mark is entered;

- 2) The database is automatically generated and the contents can't be altered;
- 3) If the test was not conducted this time, the corresponding item will be blank.


Database Information

Date/Time	Ear Mark	BF (mm)	Pregnant	Remarks
2017.05.22 10:18	4125	16.3	No	
2017.05.22 10:19	123	18.6	Yes	
2017.05.22 10:31	4568	18.2	Yes	
2017.05.23 09:23	123	17.8	Yes	
2017.06.03 17:26		18.8	Yes	
2017.06.03 17:27	982000374372458	19.1	Yes	
2017.06.15 14:52	412	17.9	Yes	
2017.06.20 17:11	4789	17.3	No	
2017.06.20 17:26	52389	18.8	Yes	

3.5 Transmission of Saved Files

The files saved by SmartScan are saved in a certain location of the smart device in four formats, images in jpg, texts in pdf, forms in excel and videos in mp4.

Search path:

All saved files can be sent or transmitted through the following functions of the smart device:

- 1) Bluetooth
- 2) WIFI
- 3) USB (the digital cable is packed in the white small carton box)

Reminder:

To transmit the saved files via WIFI, the SmartScan should not be in operating mode.

4 Backfat Measurement


The SmartScan can automatically identify and measure BF and save the measurement results in the database.

Reminder:


Backfat measurement is designed in premium version A and standard version B only. There is no the function in the basic version C.

4.1 Steps for backfat Measurement

- 1) Press the PT/BF key to switch system to BF mode. The characters BF will be display on the right upper corner of the screen in steady of PT.


- 2) Place the probe on P2 point of the animal.
Obtain a correct ultrasound image and then freeze it.

BF Measurement

- 3) Obtain a correct ultrasound image and then freeze it.
- 4) After the image is frozen, to press on the upper area of the screen gently, a green line will be appeared automatically on the image. The end of the green dotted line is the lower end of the BF line. The length of the green line represents the thickness of BF.

Meanwhile, the value of BF thickness will be immediately displayed on the left of the screen in the following format:

**Important:**

High configuration and standard configuration models can automatically identify BF on PT mode, but the data measured on BF mode is most accurate.

The BF data measured on PT mode may have a deviation of 3 to 5mm.

4.2 Save the result of backfat measurement

The measurement result of BF will be saved in the datas of the system automatically when its number appeared on the screen.

Save path: RAM directory\MYJC\report\

Save format: excel.pdf

Name: report

How to transfer the datas?

Please reference Chapter 5 below.

5 Storage & Backup of the saved files

The images, videos and reports and corresponding databases during diagnosis can be saved by the smart device. The images, videos and reports and corresponding databases saved can be connected with other smart devices via Bluetooth or WIFI and transmitted wirelessly via QQ, WeChat or other online apps.

Reminder:

To transmit the saved files via WIFI, the smart device connected to the SmartScan must be in non-operating mode.

5.1 File Saving Formats

- Image Saving Format: jpg
- Report Saving Format: pdf
- Video Saving Format: mp4
- Database Saving Format: excel

5.2 File Saving Path and Folder

All files are saved in the MYJC folder under local root directory:

File Management/Local/RAM/MYJC/

- ☐ Image
- ☐ Report
- ☐ Video

5.3 Data Backup

The data saved in the smart device should be backed up from time to time to prevent data loss.

5.4 The safest way to backup

Please to have the tablet connected with a external computer, which will to backup all the datas of SmarScan, with the digital cable, which is provided by the manufacturer and packed in the small white box in the package. (reference the Illustration below).

You are strongly to backup all th datas and have the file named as “backup13102017” .
It is better to backup one time each week.


6

Power Charging

The battery level icon of SmartScan can be found in the lower right corner of the operating interface of the smart device. When the battery is fully charged, there will be 3 bars.


Always note the change of the icon. Charge the SmartScan when there is only one bar left.

The SmartScan battery is built-in. Included with the SmartScan is a battery adapter with DC031A port.

Technical Parameters of Adapter: 4.2V, 2A, VAC100-240, 50/60Hz

The battery level when delivered from the factory is 10% only. Make sure that the battery is fully charged before first use.

At the end of the host (probe) is a transparent waterproof rubber cap. Protected by the cap is the DC-031A charging port. Open the cap and insert the charging cord. The indicator will turn red, indicating that the battery is being charged. When the indicator turns green, the battery will be fully charged and the charger can be

disconnected.

Charging Indicator Descriptions:

Color	Display	Meaning
Red	Constant	The SmartScan is being charged.
Green	Constant	The SmartScan is fully charged.

Battery Replacement:

The battery used by SmartScan is standard general battery (model NCR18650B) with protective circuit. It can be purchased on the market.

To ensure waterproof status, SmartScan adopts a built-in battery design. Please contact SONOWIN or a local dealer for professional battery replacement.

Unauthorized battery removal or replacement is prohibited. Any damage arising from unauthorized battery removal will be the responsibility of the user.

Follow the instructions provided in respect of charging and other matters of the smart device.

7 Maintenance

7.1 Routine Check and Maintenance

- Check the housing for cracks.
- Check the white oiler for cracks.
- Check the cable for cracks or other signs of damage.
- Check the battery level and see whether it's below the threshold limit value. If so, charge the battery in a timely manner.
- Check the battery charger and make sure that it's working.
- Do not drop the probe onto the floor or hard object. Do not bump the probe against other objects, the white oiler in particular.
- Do not expose the product to a dusty and humid environment to prevent degradation of product performance and reliability.

- Clean and dry the product and put it back into the original package when not in use for a long period of time.
- Properly store the product in the warehouse and comply with the storage requirements of the product.

**Warning:**

- The probe and the oiler are applicable to this product only. All consequences arising from unauthorized replacement will be borne by the user.
- The IPX7 waterproof housing of this product is one piece molded. Please contact SONOWIN for repair or replacement of any damaged part. Users should never repair it by themselves.

7.2 Cleaning and Disinfection

- The probe should be cleaned after each use and disinfected when necessary.
- With an IPX7 waterproofing rating, this product can be washed in water, but attention should be paid to the following:
 - Turn off the power.
 - Check whether the waterproof rubber cap of the battery port is properly closed.
 - Wash the product with your hands or a towel.
 - Do not clean, scrub or scratch it with hard objects or steel brushes.
- If disinfection is needed, wipe it with a cotton ball or a piece of soft cloth dipped in alcohol.
- Never clean this product with chemical liquids.

Exclusive Distributor in Europe
Schippers Europe BV
Postbus 122
5530 AC Bladel
The Netherlands
Web Site: www.schippers.eu
Chamber of Commerce: 17042016

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
Exclusive Distributor in Europe

MS Schippers
Postbus 122
5530 AC Bladel
The Netherlands
Web Site:www.schippers.eu
Chamber of Commerce: 17042016


How to set the backfat reference Line on SmartScan ?

It is known that different sows have different backfat thickness in different breeding stages. By setting the target value in advance on SmartScan before we check backfat with this wireless ultrasound, we can intuitively know whether the backfat of the pig under inspection is up to the standard during ultrasonic inspection.

SmartScan presets two different working modes: one for pregnancy Test and another is for backfat. Normally, the Pregnancy Test is performed when the machine is switched on. The scanning depth is at 160mm and PT (Pregnancy Test) character is displayed in the upper right corner of the screen. To press [PT/BF] key when the machine is in real-time, to change the system to "backfat mode". The scanning depth is at 90mm when it is in Backfat Mode, and the characters in the upper right corner are replaced by BF (Back fat).


Pregnancy Mode

Backfat Mode


Where is Reference Line?

There is a horizontal dotted line automatically appeared on the screen and the number 18.0mm displayed at its right end. This dotted line is "Reference Line" for backfat measurement.


How to move and set Reference Line ?

1) In order to prevent accidental contact, you are asked to use your finger belly to move and set this reference line.

2) Tap on the reference line with your finger belly and change it from white to yellow. Then move the dotted line up and down. The depth scale number on the right side of the dotted line will change as it moves up and down.

3) When the dotted line moves to the target value, to have your finger leaves the screen, the line color changed back to white. That means the reference line is set well.

The upper and lower limits of the reference line

5.0mm~35.00mm

Minimum accuracy of each movement

0.1mm.

Reference Value for Line:

	Very Thin	Thin	Optimal	Fat	Remarks
General Scoring System (condition from 1 to 5)	<15 mm	15-18 mm	18-20 mm	20-23 mm	