Pregtector doppler pregnancy detector

1609999

NL	Pregtector doppler blauw, schuine sonde	
	Gebruiksaanwijzing	3
EN	Pregtector doppler pregnancy detector	
	Manual.....	8
DE	Pregtector Doppler, abgewinkelte Sonde	
	Gebrauchsanweisung.....	13
FR	Echographe Pregtector Doppler bleu	
	Mode d'emploi.....	18

Pregtector doppler blauw, schuine sonde

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Pregtector doppler pregnancy detector

1. Checking the instrument	10
2. Before using the PregTector.....	10
3. Using the PregTector	10
A) First attempts on Pregnant sows.....	10
B) Early Pregnancy Diagnosis.....	11
C) Negative Diagnosis	12
4. A suggested pregnancy detecting routine.....	13
5. Pregnancy detection before the first return at 21 days	13
6. Using the doppler pregtector during oesterus	14
7. Parvo virus and s.m.e.d.i.....	14
8. Using the dopller pregtector during farrowing	15
9. Pregnancy detection of outdoor sows using the doppler PregTector	15
10. Battery changing	17
11. Guarantee.....	17

1. Checking the instrument

Unpack the PregTector carefully, keeping the packing materials safety in case they should be needed at any time in the future for returning the instrument for service. In the outer packaging shows any sign of damage in transit, please immediately inform or the dealer from whom the PregTector was purchased.

The following are the standard parts which should be in the package:

1. Instrument carry case
2. Headset
3. Probe containing circuit and battery
4. Plastic bottle of coupling medium (liquid paraffin)
5. Instruction manual
6. Instruction cassette

In the unlikely event that any parts are missing or damaged, please immediately inform or the dealer from whom the instrument was purchased.

2. Before using the PregTector

Plug the headset into the lead coming from the probe, put on the headset and adjust it to be a comfortable fit. Turn the volume controls to maximum. Release the switch and put a quantity of coupling medium on the probe-face. Apply the probe to your own chest, next to the skin, over the heart, which will normally be between the left nipple and the breastbone. Depress the switch, and your own heartbeat should be clearly heard. Should the "click" of the instrument being switched on not be heard, not the low hiss indicating that it is switched on, change the battery. It is possible for batteries to be damaged during transit. Always use a alkaline battery.

Before trying to use the PregTector on a sow, listen to the instruction cassette, and read this manual.


3. Using the PregTector

A) First attempts on Pregnant sows.

Don't begin by trying to test a doubtful gilt just four weeks served, but start if possible in the farrowing house, on a quest sow waiting to farrow. This will be easier because the uterus is of course quite big by now, and also because there will usually be more space in the farrowing house than in sow stalls.

Always try to test animals when they are as quiet as possible. This will minimize noise caused by the sow's movements, which will make identification of the sounds more difficult, especially for an inexperienced operator. Try also to avoid as far as possible the first four hours after feeding, to minimize noise caused by digestive process. There is no need for the sow to be standing up, as with some other detectors. If there is no sow waiting to farrow which is quiet, it is best to wait until there is. It is rarely if ever essential to test any animal at any particular time and as with many other tasks you will do with animals, pregnancy testing is much on non-stressed, calm sows.

Place the headset comfortable on the top of the head and cover the probe-face with coupling medium. Place the probe on the sow's flank, in front of the leg, on the belly, approximately in line with the second teat from the rear. (see illustration). Direct the probe upwards, imagining that the beam of ultrasound will pass through the centre point of the sow and emerge on the opposite side of her backbone. The uterus runs roughly from front to back, so a small vertical 'scanning' movement will usually pass the sound beam over the uterus. On this heavily pregnant sow, you probably hear the fast pulse of the umbilical cord, or even a fetal heart, immediately the unit is switched on. Adjust the volume control for maximum comfort and clarity. Listen to these sounds recorded on the instruction cassette in order to be sure of recognizing them. Having found one of them, search for the other by changing the angle of the probe. Apply firm but not excessive pressure to the probe. One result of pressing too hard against the sow is that she will be disturbed by this and move around more, causing unnecessary noise in the headphones. In case of a sudden movement or other unexpectedly loud noise, immediately release the on/off switch to silence the detector.

Be careful not to move the face of the probe across the sow's skin. This will cause noise in the headphones and also remove the coupling medium from the probe. You can move the probe and loose belly skin together and cover quite a large area of the abdomen without causing any noise or losing contact. Should the probe become removed from the sow, either because you want to try listening in a different site, or because movement of the sow has necessitated removal, you must always re-apply coupling medium before replacing the probe. Coupling medium has the effect of removing air from between the probe and the skin. Sound is absorbed by air, and the 'volume' of the ultrasound being used here is so small that any air it passes through will greatly affect the performance of the instrument.

It may be difficult to find the uterine artery in this heavily pregnant sow, because it runs on the upper surface of the uterus, which by now is swollen with its contents, and it could be beyond the range of the detector. You could however try for a minute or two to find it. It will be clearly different from the fetal pulses you have already heard because the sow's heart rate will be only about half as fast as that of the fetuses. If you cannot find the sow's pulse, don't worry, as the uterine artery is largely meaningless at this advanced stage of pregnancy, as long as you can hear one or more fetal pulses. Carry on familiarizing yourself with the instrument for a while before moving to the dry sow house to check a sow around the mid-point of gestation.

B) Early Pregnancy Diagnosis

Although you can test sows from as early as 14 days after service or insemination, your first attempts should be with sows around 6 weeks pregnant. Again, choose a lying sow to begin with as this will make testing easier. Cover the probe with coupling medium and place it in the apex of the triangle of belly skin, as illustrated. Don't turn on the PregTector until you are ready to listen. When you do, listen for the whip-like, whistling, pulsing sound of the uterine artery. This will be much slower than the fetal pulses you were listening to in the farrowing house, as is shown on the instruction tape. Once you are sure that this sound is present, the diagnosis has been made. Now move on to other sows, gradually testing earlier and earlier in gestation. You will quickly learn to make a positive diagnosis, as the sound of the uterine artery is so distinctive. As you gain in confidence, you will find that you only need to hear two or three pulses from each sow, and testing will be very quick.

C) Negative Diagnosis

In the nature of the way the Doppler PregTector Works, where you are making a positive diagnosis by recognising a positive signal, negative diagnosis is a little more difficult. If the sow is not pregnant, the characteristic sounds you have been listening to are absent. You will clearly have to spend a little longer making a certain as you that what you are listening for is definitely not there, and that your lack of success is not due just to looking in the wrong place. There are some rules to follow when confirming a sow open or not-in-pig.

1. Don't make a negative diagnosis on just one test. Before confirming a sow is not pregnant, test her at least twice, with an interval between tests of as long as can be sensibly arranged, ideally a week.

2. It is important to know the service date when making a negative diagnosis. It is obviously quite easy to confirm a negative if the sow should be ready to farrow, because at that time it should be very easy to find one or other of the vital sounds of pregnancy, as the uterus presents a large and easy-to-find target. But the uterus in early pregnancy is small, and can be hidden behind other organs, which can block the passage of ultrasound, or else place the uterus beyond the range of the transmitted ultrasound.

3. If the sow is lying down, make her stand up. It can happen that sleeping sows are disturbed when you begin testing them, and so sometimes you test sows that have just stood up. In early pregnancy the uterus is small and can be displaced behind other organs, for example, the intestine. When the sow stands up, the various structures in her belly do not immediately fall into their normal positions, and it can happen that you don't hear the uterine artery because the uterus is temporarily not in its usual place.

If you test that sow again five minutes later you will quite possibly hear the uterine artery loud and clear the moment you switch on the PregTector, and so you can confidently mark her as pregnant.

4. To erroneously make a negative diagnosis is not a serious mistake, unless you immediately cull the sow as a result. Your attention will have been drawn to that sow, and she can be retested at any time thereafter to confirm that she is definitely not in-pig. In the worst case, you would certainly keep the sow for at least three weeks to check her for return to oestrus and test her again before culling. If she was pregnant all the time, she will now be three weeks more advanced in gestation, and the uterus correspondingly larger and therefore more difficult to miss. At least your attention will have been drawn to that animal as being a potential problem.

On the other hand, a false positive diagnosis, where you confirm a sow pregnant when she is not, is the most expensive error you can make. There will be a tendency to ignore that sow, or at least not pay any special attention to her, and it is not unknown for such sows to be kept for the whole gestation period, and not be discovered non-pregnant until they fail to farrow. When the cost of labour, finance, feed and loss of profit on the lost litter is calculated, you will quickly understand the value of accurate pregnancy detection.

5. Making a negative diagnosis becomes much easier with experience. Once you know that you can find positive pregnancy, you will be confident that the lack of the sounds can definitely be taken to mean that the sow is not-pregnant. A practise which does not help negative diagnosis is just to pregnancy test the doubtful sows. Generally, the mere fact that they are doubtful means that they will be more difficult. Test all sows systematically, so that you become confident in using the PregTector on the easy ones as well. Most of them are easy, so the fact that you occasionally come across one which is not is straight away a possible indication of a problem.

4. A suggested pregnancy detecting routine

A pregnancy detector will not take the place of traditional animal husbandry, but it can be used in conjunction with it to improve efficiency in the farm. Management decisions in pig breeding, as in any other venture, will be better in proportion to the accuracy, validity and up-to-dateness of the information needed to make them.

Sows should ideally be kept in small groups in the service area for 4-5 weeks after service or insemination. They should be visually checked for return to oestrus at 21 days, and they should be introduced to the boar as an extra insurance to pick up returns.

The Doppler PregTector can be used from 14-15 days, but a more usual time is at about 4-5 weeks.

The suggested routine is as follows:

1. Check sows systematically and regularly on a routine basis each week. Choose a day when there is less activity on the farm. (Not the day when pigs are weighed, or loaded, or when sows are weaned). Do the pregnancy testing at a quiet time of day, ideally not within four hours after feeding, and certainly not just before feeding them, when the sows are usually at the noisiest. Stick to this routine rigidly.
2. Test sows in their fifth week, that is, they will be between 28 and 35 days. Sows which are confirmed pregnant are no problem for the time being, but those where nothing is heard must be tested again later. Check the exact service date. The important thing is to get them back to the boar, or re-inseminate at the next return, at 42 weeks or so. Therefore if the sow is 28 days served, you can just make a note to test her again next week. But if she is at the end of her fifth week, she will need to be tested individually two or three days later so that if a negative result is obtained again, she can be prepared to be served or inseminated again at the appropriate time.
3. Some farms will routinely test a second time during the week after the second return, as it is possible for a sow to reabsorb the pregnancy up to about 45 days, and the most likely times are at the first and the second return days. Others will leave the second test until about 10 weeks, which is more or less the mid-point of the pregnancy. If testing later than 10 weeks you should listen not only for the uterine artery but also for the faster pulse of a fetus.
4. A Major advantage of the Doppler PregTector is that pregnancy testing can be carried out at any time after 25 days. Use the instrument at any time you have doubts about a particular sow.

5. Pregnancy detection before the first return at 21 days

The physical changes which cause the uterine artery to become enlarged in early pregnancy are first detectable in most sows from day 14 to 15 after service with the PregTector. Trials at the University of Göttingen in Germany (Holz et al Top Agrar Jan 1994) have found that positive diagnosis was possible in more than 90% of sows from the 14th day. Because of the uterus is very small at this time, such a diagnosis will require more experience than testing later in pregnancy.

It must also be understood that any diagnosis made before the first return to oestrus at 21 day must be reconfirmed later, as there is a significant possibility of the sow losing the pregnancy,

probably by re-absorption, during that first cycle. However, such a test can never-the-less be very useful in some circumstance

This could be that there is a shortage of accommodation in the service area, and the breeder would for the sake of convenience like to move sows to the dry sow accommodation rather sooner than might be strictly recommended. The ability to be sure that a sow has conceived, and therefore will at least not return at 21 days (even though she may do so later, perhaps at 6 weeks) will ensure that only pregnant sows are moved.

Another reason could be to try to identify the root causes of a temporary infertility problem; if sows are conceiving but then returning to oestrus later, then the problem has been isolated to the sow, but if sows are not becoming pregnant at all, this will indicate that the fault could lie with the boar, or there could be some problem of artificial insemination.

It is unlikely that most breeders will want to test sows as a matter of routine as early as 14 days, but the ability to do so when necessary is none the less useful.

6. Using the doppler pregector during oestrus

Occasionally during oestrus, blood-flow in the uterine artery will be detectable. This is more likely in gilts and small sows. The characteristic whistling sound of true pregnancy will be absent, and the artery will have a lower-pitched, quieter sound, and there will be a period of silence between the pulses. The operator should be aware of this, but in practise it will not be a disadvantage. If there is any suspicion that the sow is on heat rather than pregnant, then other checks, such as visual examination, trying the standing reflex and introducing the sow to a boar should be put in hand immediately. If this is not possible, then that sow should be rechecked with the PregTector about 4 days later. If she was on heat, the uterine artery will have ceased, and if pregnant, will have become stronger.

7. Parvo virus and s.m.e.d.i.

Whilst the fetal pulse in the umbilical cord may be heard as early as 31 days. And the fetal heart from 40 days, their absence is not really significant until the second half of pregnancy. A sow with parvo-virus or S.M.E.D.I. can be carrying dead fetuses. She is, in reality, pregnant, but not productively so, which is what the breeder needs to know. The uterine artery can sometimes be heard in such sows, although it will not be as easy to find or as clear as it would be in a normal pregnancy, but since both these viruses tend to strike well into the second half of pregnancy, it is the absence of fetal sounds which must be the means of diagnosis. A useful aid to diagnosis of either of these problems is to compare doubtful sows with others of the same or similar gestational age, as the difference is usually quite apparent.

8. Using the doppler pregtector during farrowing

The ability to confirm the presence of live fetuses in the uterus can be very useful in a difficult farrowing. Whilst it is nearly always possible to manually examine the uterus, this should be avoided if at all possible as, apart from being unpleasant both the breeder and the sow, it can lead to uterine infection. There are times when it cannot be avoided, but the Doppler Pregtector can help to avoid its use for exploratory purposes.

Occasionally a sow will deliver a small number of pigs, followed by some cleansing, and uterine contractions will stop. The placenta, which has been nourishing the fetuses throughout pregnancy, will begin to break down at the moment the first fetus is expelled.

Subsequent pigs must be born within 6 to 12 hours of the first in order to survive, and in these circumstances, if nothing is done, any further pigs will be stillborn within the next 12 to 24 hours. If this is suspected, it will take only a minute or two to check the length of the sow's belly with the Pregtector. The uterine artery will be much in evidence, but this is no longer significant, as it will continue to be detectable for up to three weeks after farrowing, as the uterus returns normal following the pregnancy. However, if the faster fetal pulse can be heard, this is definite proof that a live fetus is in the uterus. Steps can now be taken, for example manual intervention, or an injection of oxytocin to restart contractions, to deliver the remaining fetus or fetuses as quickly as possible.

9. Pregnancy detection of outdoor sows using the doppler Pregtector

Pregnancy detection of outdoor sows is no less important than their indoor counterparts. As with most procedures carried out on outdoor sows, it can be a little more difficult, but steps can be taken to smooth the way.

Doppler pregnancy detectors like the Pregtector work by detecting the vital movements of pregnancy: blood flow in the uterine artery and umbilical cord, and the movements of the fetal heart itself. These movements cause small changes in the frequency of returned echoes, which are amplified by the Pregtector and which closely resemble how we might expect these structures to sound. Movements of the sow herself are much greater than the movements associated with pregnancy, and will inevitably tend to 'swamp' what we are looking for.

The Pregtector has a system of audio filters which 'target' the instrument at the range of frequencies we are looking for. This results in a definite tendency to find the sounds more quickly, which is very helpful with unrestrained sows. But some of the frequencies associated with sow movement are within the targeted range, and a part of the noise results from the changes in pressure between the probe and the sow caused by movement of both sow and operator, and these cannot be filtered out without severely limiting sensitivity of the detector. Therefore it remains important that the sow should be as quiet as can be arranged.

Many outdoor breeders are successfully using the Pregtector, and their experience can be useful to others starting out on this road.

First, one or two main points to mention:

The uterus becomes bigger, and therefore easier to find, day by day. The later you can leave diagnosis the easier it will be to locate. Whilst the increased blood flow in the uterine artery can be

heard as early as 14 days, unless you really need to know as early as that, leave it until later. The important date by which to know is the second return at 42 days, so leave testing until the last week before this if you can.

Sows can be tested on either side, either standing or lying (sitting is sometimes difficult).

Of course, each farm is different, so each will find slightly different ways of optimising performance of the PregTector. Some farms find that testing can be done in the paddocks, just by approaching the sows, if you are able to do this while they are lying, testing will be as easy as testing indoor sows.

The key to simple testing is find the right time in the sows' day when they are most relaxed. This generally will be two hours or more after feeding. Continue sows to be tested to a small paddock or central handling area, and allow them to settle. Ideally, they will be in groups of not more than 6 or so, at least for the purpose of testing, even if they are living in larger groups. If you don't have an area where sows can be restricted, at least to a degree, it may be necessary to use a trailer, or rig up some kind of moveable pen perhaps using hurdles.

If sows will stay in one place, without being persuaded to do so by a bucketful of nuts, this is generally better. The action of eating causes quite a lot of movement of abdominal muscles, which can be noisy. However, this is quite likely to be impossible. You will find they settle more quickly if they have been fed.

Just as important probably, is to test systematically, so that unnecessary testing is avoided, and to test regularly. It is far better for the operator to test a little and more often than to try to have a marathon testing session every once in a while. Familiarity will quickly bring confidence in the PregTector, and testing will be a routine matter. It is also a fact that sows become used to being tested and hinder the operator less.

As far as the use of the PregTector itself is concerned, there is no difference in technique from that described in the Operators Handbook for restrained sows. Be aware that excessive pressure of the probe on the sow will cause her to move more, and that mud on her flanks may have to be removed.

It is rarely, if ever essential for a particular sow to be tested at a particular moment. So if the going is difficult, leave it for a while until conditions improve.

Finally, keep a sense of humour, you may need it!

10. Battery changing

The PregTector has a low battery warning in the form of a beep tone. This will start to sound once the battery voltage has dropped to a predetermined level, and the speed of the beep will increase as the voltage drops.

The PregTector will remain usable for a time after the tone sounds, but you should ensure that a new battery is fitted at the earliest opportunity. The battery is located in the probe at the opposite end to the probe face. To change the battery remove the end-cap by twisting in out, and shake out the battery with its connector. Always replace with a new battery of the same type. Use only alkaline batteries as these will give far better performance. Battery life will depend on the volume level chosen. In case of any suspected problem with the detector, always change the battery before considering returning the instrument for service. Rechargeable ni-cad batteries can also be used, but unless you have difficulty in obtaining new alkaline batteries, or already have suitable rechargeable cells, it is recommended that alkaline batteries are used, as their performance will be better. Batteries should ideally be bought from a store with a reasonable turnover of stock, because the shelf-life of alkaline batteries is limited to about 2 years.

11. Guarantee

The Doppler PregTector is guaranteed for a period of 3 years from date of invoice, against defective parts and manufacture. Instruments will be repaired free of charge provide that:

- a) The unit is returned, carriage paid to the dealer from whom it was purchased.
- b) The unit had been damaged by misuse, by maltreatment, including damage caused by animals, or by unauthorized attempts at repair.
- c) Batteries are specifically excluded from the warranty.

12. Trouble shooting

You may be able to save yourself time and trouble by carrying out simple checks and repairs yourself. This is particularly the case where the PregTector is completely 'dead'. There are 2 likely causes of this.

1. Check that the battery connector is making good contact with the two terminals on the battery. Check especially that the female (+ve) cup on the connector side is properly gripping the male pole on the battery. If not, a gentle squeeze with pliers may solve the problem.
2. Make sure that the volume controls on the headset are full open. You may be able to check that the connectors between the headset and the probe are functioning correctly if you have another headset on a piece of audio equipment which had the same 6mm (1/4) jack plug. If the 'hiss' when the on button is depressed can be heard with the test headset, but not with the original, the problem is probably in the headset. But do not attempt to use the PregTector, except for this short test, with any headset without first ensuring that it has an impedance value of 8 ohms.

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**Schippers Italia Srl**
**Schippers Schweiz GmbH**
**Schippers UK Ltd.**
**LLC Schippers Russia**
**Schippers Canada Ltd**
7102 52nd Street Bay # 18 • Lacombe, AB • T4L 1Y9