

CONTRAfluran™

Narkosegasfilter

Der CONTRAfluran - MS Narkosegasfilter adsorbiert volatile Anästhetika. Bei fachgerechtem Einsatz schützt er den Patienten und das Personal im Behandlungsraum.

Bitte beachten:

1. Gebrauchsanleitung lesen.
2. Nur einen Narkosegastyp je filter einsetzen
(Gastyp bitte auf dem Etikett vermerken)
3. Lagertemperatur: -5 °C bis + 35 °C
4. Lagerung in trockenen, staubfreien Räumen
Relative Luftfeuchte: max. 70%

MS Schippers

www.schippers.eu

CONTRAfluran™

Narkosegasfilter

Schippers GmbH

Kölner Straße 62
47647 Kerken (D)
Tel: 02833-92360



Artikelnummer: 3200136
Medizinprodukt: Klasse 1
Made by SeoSys Medical GmbH

MS Schippers

www.schippers.eu

CONTRAfluran™

Narkosegasfilter

The CONTRAfluran Anaesthetic Gas filter adsorbs volatile anaesthetics. Its use protects staff and patients.

Please take note of:

1. Read the manual before use.
2. Use only one type of anaesthetic gas per filter.
(select the type of gas at the cap label)
3. Storage temperature: -5°C to + 35°C
4. Storage rel. humidity: max. 70%

MS Schippers

www.schippers.eu